



A.T.S.[®] Contour Cuffs



Versatility and flexibility for both surgeon and patient



No two people are the same.

Each patient is different, which makes each procedure unique. That's why Zimmer offers the versatility and flexibility of A.T.S.® Adult Contour Tourniquet Cuffs.

A.T.S. Contour Cuffs are designed to fit most adult patients for arm, thigh, lower leg and Bier Block (IVRA) applications. In addition, Contour Cuffs use unique, patented pivoting straps with a contour design that adjust to a wide range of limb shapes and sizes. The shape is optimized for conically-shaped limbs.

Because they are wider than other cuffs, A.T.S. Contour Cuffs may provide improved application of pressure to the operative limb.¹ And because one cuff fits all patients for a particular limb, Contour Cuffs offer a simple solution for inventory control.

Depend on A.T.S. Contour Cuffs to meet any procedural demand without sacrificing patient comfort and safety.





SUPERIOR FIT AND PERFORMANCE

A.T.S. Contour Cuffs have a slim, low-profile design which provides superior fit and performance on tapered limbs.¹ The cuffs are durable and easily cleaned for multiple uses. Three Contour Cuff styles are available.



CONTOUR THIGH CUFF

Offered in both regular and large sizes to fit thighs from 9 to 39 inches (23 to 99cm) in circumference. Available in single port, single bladder and dual port, single bladder options.



CONTOUR UPPER LIMB CUFF

Fits upper limbs from 9 to 21 inches (23 to 53cm) in circumference. Available in single port, dual port, single bladder and dual bladder options.



CONTOUR LOWER LEG CUFF

Fits lower legs from 9 to 17 inches (23 to 43cm) in circumference. Available in single port, single bladder and dual port, single bladder options.

The Zimmer comprehensive A.T.S. Contour Cuff line provides multiple options to help the clinician meet each patient's specific needs.

1. Younger ASE, McEwen JA, Inkpen K. Wide contoured thigh cuffs and automated limb occlusion measurement allow lower tourniquet pressures. *Clin Orthop*. November 2004;(428):286-293

A.T.S. Contour Cuffs

- Reusable
- Easy to apply
- Low profile
- 90 degree angle ports
- Latex free, durable
- Various lengths - 9 to 39 inches (23 to 99cm)
- Available in Dual Port, Single Port, Dual Bladder and Single Bladder designs

Contour Cuffs*

The cuffs listed below include one sleeve unless otherwise indicated.

Dual Port, Single Bladder

Prod. No.	Description	Limb Size
60-7450-002-00	Contour Upper Limb Cuff	9 to 21 in.
60-7450-003-00	Contour Thigh Cuff	12 to 30 in.
60-7450-004-00	Contour Lower Leg Cuff	9 to 17 in.
60-7500-008-00**	Large Contour Thigh Cuff	30 to 39 in.

Dual Port, Dual Bladder

Prod. No.	Description	Limb Size
60-7450-011-00	Contour Upper Limb IVRA Cuff	9 to 12 in.

Single Port, Single Bladder

Prod. No.	Description	Limb Size
60-7350-002-00	Contour Upper Limb Cuff	9 to 21 in.
60-7350-003-00	Contour Thigh Cuff	12 to 30 in.
60-7350-004-00	Contour Lower Leg Cuff	9 to 17 in.
60-7600-008-00**	Large Contour Thigh Cuff	30 to 39 in.

Single Port, Dual Bladder

Prod. No.	Description	Limb Size
60-7350-011-00	Contour Upper Limb IVRA Cuff	9 to 12 in.

*Sold Non-Sterile

**A.T.S. Contour Cuffs are 5.75 inches wide, except 60-7500-008-00 and 60-7600-008-00, which are 4.75 inches wide and do not come with a sleeve.

Limb Protection Sleeves

- Specific designs for Contour Cuffs
- Shown to be more effective at eliminating wrinkles and pinches in the skin surface²
- Color coded to match the accompanying cuff
- Does not clog cuff closure material

Limb Protection Sleeves*

For Contour Cuffs (Box 10)

Prod. No.	Description
60-7920-001-00	Contour Upper Limb Cuff Sleeve
60-7920-002-00	Contour Thigh Cuff Sleeve
60-7920-003-00	Contour Lower Leg Cuff Sleeve

*Sold Non-Sterile

Esmark Bandages

Latex Free Sterile (Box 20)

Prod. No.	Description
60-2694-011-00	3 in. x 9 ft. (8cm x 3m)
60-2694-012-00	4 in. x 9 ft. (10cm x 3m)
60-2694-013-00	5.5 in. x 9 ft. (14cm x 3m)
60-2694-014-00	3 in. x 12 ft. (8cm x 4m)
60-2694-015-00	4 in. x 12 ft. (10cm x 4m)
60-2694-016-00	5.5 in. x 12 ft. (14cm x 4m)

Sterile (Box 20)

Prod. No.	Description
60-2694-001-00	3 in. x 9 ft. (8cm x 3m)
60-2694-002-00	4 in. x 9 ft. (10cm x 3m)
60-2694-003-00	6 in. x 9 ft. (15cm x 3m)
60-2694-004-00	3 in. x 12 ft. (8cm x 4m)
60-2694-005-00	4 in. x 12 ft. (10cm x 4m)
60-2694-006-00	6 in. x 12 ft. (15cm x 4m)

2. McEwen JA, Inkpen K. Tourniquet safety preventing skin injuries. *The Surgical Technologist*. August 2002;7-14.

Contact your Zimmer representative or visit us at www.zimmer.com

