



World Leader in
Orthopaedic Products

Zimmer Orthopaedic
Braces and Supports





About Zimmer

Justin O. Zimmer started Zimmer Manufacturing Co. in 1927 with a line of aluminum splints which were introduced at the American Medical Association meeting held that year in Washington, DC. Since that time, Zimmer has been a leader in the orthopaedic industry. Today, Zimmer, Inc. employs over 6,000 people worldwide and maintains a leadership position in virtually every segment in which it competes. Our market leadership in orthopaedics and reputation for innovation, superior service and excellent product quality has been built over the past 75 years. Our interest in orthopaedic braces and supports originated from these roots.

We have a worldwide network of distributors and representatives eager to help you and your patients. In an independent survey of U.S. orthopaedic surgeons, Zimmer was ranked as the number one company in 10 of 11 categories related to products and service. The company's sales force was once again named as most responsive and most knowledgeable. The Sales Associates undergo extensive training to help ensure that they have a thorough knowledge about orthopaedic products, procedures and treatments providing them a better understanding of orthopaedics and the needs of orthopaedic institutions. There are benefits associated with working with a recognized sales team and the excellent services and resources provided through Zimmer Orthopaedic Surgical Products. Through unparalleled service, exceptional products and education, we put "*Confidence In Your Hands.*"

At Zimmer, we consider the quality of our fabric and closure materials as important as meticulous design and construction. After selecting the highest-quality materials best suited to each product application, we maintain strict adherence to fabric specifications through periodic testing, to further ensure product performance and patient comfort. No stiffening agents are ever added to Zimmer product to create the appearance of extra body and durability. Our materials are selected to retain their original characteristics through many washings. Each and every product is delivered to you with Zimmer pride.

The quality inherent in our manufacturing process and our committed staff enables us to offer our customers exceptional products at reasonable prices. While we understand that cost containment pressures are growing and that product pricing has become an important issue, you can still provide your patients with the best products available within affordable limits by selecting a *Zimmer*[®] Orthopaedic Brace or Support. Your Zimmer representative will strive to provide you with competitive pricing for our high quality braces and supports.

Special Needs, Custom Orders

We offer a comprehensive, full-line selection of quality orthopaedic braces and supports including pediatric sizes. The diversity and flexibility of the product styles as well as availability of a large range of sizes allows you the opportunity to select the supports most appropriate to your patients' specific needs.

Because patients come in all sizes and shapes, we understand that they may sometimes have special requirements. We will try to accommodate your special requests, so that each patient has the best fit possible. Through a special ordering process many unique sizes that address specific patient deformities or unusual anatomies are available. **Extra-large and extra-small sizes are often also available through this special order process. Some of these are listed in our catalog, so that you can reference a specific catalog number. These items are denoted by a ♠. We also offer customer-specific logos on a limited list of products. Many of these are noted in our catalog. Please note that custom products cannot be returned for credit.** For more information contact Zimmer Orthopaedic Surgical Products at 800-321-5533.

Reimbursement

Zimmer does not warrant the appropriateness of reimbursement codes. The provider is solely responsible for billing products under the correct reimbursement code.

Ordering

Please contact your local Zimmer representative or call Zimmer Customer Service at 1-800-348-2759.

Visit Our Website

www.zimmer.com



TABLE OF CONTENTS

Cervical Collars	6
Upper Extremity/Shoulder/ Clavicle Supports	10
Wrist Supports	17
Abdominal/Back/Rib Supports	20
Knee Immobilizers	23
Knee Braces and Supports	26
Lower Leg/Ankle Supports	31
Walkers	34
Cast Shoes/Post-Op Shoes	36
Miscellaneous Supports	39
Orthopaedic Wraps	40
Cold Therapy	41
Pressure Controlled Stockings	42
Splints	43
Casting Products	48
Pads and Positioners	52
Traction Boots	57
Pelvic Belts	58
Head Halters	59

CERVICAL COLLARS

PHILADELPHIA® TWO-PIECE CERVICAL COLLAR WITH TRACHEAL OPENING

- Anterior and posterior rigid plastic reinforcement supports limit movement
- Large tracheal opening for easy access to anterior neck region
- 2-piece design with three collar heights allow customized fit and prevention of flexion and hyperextension while maintaining neutral alignment during patient recovery
- Latex-free, non-toxic and hypoallergenic foam reduces skin irritation and provides patient comfort and compliance
- Adjustable bilateral contact closures ensure comfortable fit while allowing easy application
- X-ray, CT and MRI lucent
- Washable with mild soap and clean water
- Conforms to shoulder and mandible contours for superior stability and comfort
- Water-resistant materials allow patient bathing and aqua therapy without removing collar



ADULT			
Prod. No.	Size	Height	Neck Circumference
00-1741-005-00	Small	3 1/4 in. (8cm)	10 in. - 13 in. (25cm - 33cm)
00-1741-010-00	Small	4 1/4 in. (11cm)	10 in. - 13 in. (25cm - 33cm)
00-1741-015-00	Small	5 1/4 in. (13cm)	10 in. - 13 in. (25cm - 33cm)
00-1741-006-00	Medium	3 1/4 in. (8cm)	13 in. - 16 in. (33cm - 41cm)
00-1741-011-00	Medium	4 1/4 in. (11cm)	13 in. - 16 in. (33cm - 41cm)
00-1741-016-00	Medium	5 1/4 in. (13cm)	13 in. - 16 in. (33cm - 41cm)
00-1741-007-00	Large	3 1/4 in. (8cm)	16 in. - 19 in. (41cm - 48cm)
00-1741-012-00	Large	4 1/4 in. (11cm)	16 in. - 19 in. (41cm - 48cm)
00-1741-017-00	Large	5 1/4 in. (13cm)	16 in. - 19 in. (41cm - 48cm)
PEDIATRIC			
Prod. No.	Size	Height	Neck Circumference
00-8601-968-50	Ped	2 1/2 in. (6cm)	10 in. - 12 in. (25cm - 30cm)

ZIMMER® COOLBRACE®* CERVICAL ORTHOSIS

- A small battery-operated fan reduces heat buildup, perspiration and skin rashes associated with prolonged collar use for increased patient comfort and compliance
- Soft, washable liners for improved patient hygiene
- 2-piece design with plastic shell to maintain cervical stability and support
- Easy application
- Batteries are not included



Prod. No.	Size	Height	Neck Circumference
00-1751-005-00	Small	4 1/4 in. (11cm)	10 in. - 13 in. (25cm - 33cm)
00-1751-006-00	Medium	4 1/4 in. (11cm)	13 in. - 16 in. (33cm - 41cm)
00-1751-007-00	Large	4 1/4 in. (11cm)	16 in. - 19 in. (41cm - 48cm)

*Trademark of Leonard G. Davis and Dr. Lytton Williams



CERVICAL COLLARS

ZIMMER® PADDED CERVICAL COLLAR

- Removable pads wick away moisture to ensure better patient hygiene
- Restricts cervical spine flexion, extension, and rotation to promote patient recovery
- Injection molded polyethylene construction helps improve rigidity and durability
- Large tracheal opening for easy access to anterior neck region
- Conforms to shoulder and mandible contours to help increase comfort
- X-ray, CT, and MRI lucent
- Latex-free



Prod. No.	Size	Height	Neck Circumference
00-1741-301-00	Stout	1 3/4 in. (4cm)	11 in. - 24 in. (28cm - 61cm)
00-1741-302-00	Medium	2 1/2 in. (6cm)	11 in. - 24 in. (28cm - 61cm)
00-1741-303-00	Standard	2 in. (5cm)	11 in. - 24 in. (28cm - 61cm)
00-1741-304-00	Tall	3 in. (8cm)	11 in. - 24 in. (28cm - 61cm)

Prod. No.	Description
00-1741-309-00	Replacement Pad Set

EMS EXTRICATION COLLAR

- Clinically proven to provide effective support against flexion/extension
- Large tracheal opening for easy access to anterior neck region
- X-ray transparent
- Contoured sides place no pressure on ears
- 1-piece design stores flat for convenient storage in limited space
- Collar designed to support weight of up to 100 pounds
- Color-coded closure straps and size indicators for quick, easy and accurate selection



Prod. No.	Size	Height	Neck Circumference	Color Code
00-1676-001-00	Ped	3 1/2 in. (9cm)	17 3/4 in. (45cm)	Red
00-1676-002-00	X-Short	2 3/4 in. (7cm)	22 1/2 in. (57cm)	Green
00-1676-003-00	Short	3 1/2 in. (9cm)	22 1/2 in. (57cm)	Yellow
00-1676-004-00	Regular	4 1/2 in. (11cm)	22 1/2 in. (57cm)	Orange
00-1676-005-00	Tall	5 1/2 in. (14cm)	22 1/2 in. (57cm)	Blue
00-1676-006-00	Ped Short	2 1/2 in. (6cm)	17 in. (43cm)	Lt. Blue

PHILLY®* ONE-PIECE® EXTRICATION COLLAR

- Large tracheal opening for easy access to anterior neck region
- Restricts flexion, extension and rotation to immobilize the cervical spine
- Latex-free, non-toxic and hypoallergenic foam reduces skin irritation
- External sizing guide for easy size selection with color-coded closure straps and labels to ensure easy recognition of selected size
- 1-piece design stores flat for convenient storage in limited space
- Rear ventilation slots for fluid drainage and airflow
- X-ray, CT and MRI lucent



Prod. No.	Size	Height	Neck Circumference	Color Code
00-1750-001-00	Infant	5.7 in. (14cm)	17.66 in. (45cm)	Pink
00-1750-002-00	Ped	5.75 in. (15cm)	18.22 in. (46cm)	Aqua
00-1750-003-00	Low	6.44 in. (16cm)	22.50 in. (57cm)	Purple
00-1750-004-00	Short	7.05 in. (18cm)	22.50 in. (57cm)	Blue
00-1750-005-00	Regular	7.69 in. (20cm)	22.50 in. (57cm)	Orange
00-1750-006-00	Tall	8.31 in. (21cm)	22.50 in. (57cm)	Green

*Trademark of Philadelphia Cervical Collar Company

CERVICAL COLLARS

PHILLY®* PATRIOT®* ADJUSTABLE EXTRICATION COLLAR

- Delivers safe motion restriction without hyperextension
- Rigid 1-piece injection molded polypropylene construction for immobilization
- Oversized tracheal opening for easier access to anterior region of neck for emergency airway procedures
- Easy-to-use instructions: measure the patient, select the size, adjust height, and lock
- Essentially flat design for easy, convenient storage where limited space is available
- Radiolucent, CT and MRI compatible
- Chin support with foam padding
- Nasal cannula holders
- Open rear panel for ventilation and drainage
- Multi-height adjustment reduces the need to stock multiple sizes and allows for optimal size selection



PEDIATRIC			
Prod. No.	Size	Height Range	Neck Circumference
00-1742-001-00	Pediatric	3/4 in. - 1 3/4 in. (2cm - 4cm)	8 in. - 18 in. (20cm - 46cm)
ADULT			
Prod. No.	Size	Height Range	Neck Circumference
00-1742-002-00	Adult	1 3/4 in. - 3 1/2 in. (4cm - 9cm)	11 in. - 23 in. (28cm - 58cm)

*Trademark(s) of Philadelphia Cervical Collar Company

Pediatric version adjusts to 3 sizes

Adult version adjusts to 4 sizes

ZIMMER® SERPENTINE FOAM COLLAR

- Contoured design can be applied for slight flexion or reversed for slight extension
- Sold individually
- Choice of 2 densities – firm or medium
- Designed to provide immobilization and gentle support of the cervical spine
- Washable
- Convenient contact closures for easy adjustment and application
- Wide selection of collar heights for a proper fit
- Collars packaged with extra washable cotton blend stockinette and extension for patient convenience and comfort
- Indicated for cervical facet syndrome, muscle sprains and strains, degenerative disc disease with and without radiculopathy, discogenic pain



FIRM DENSITY			
Prod. No.	Size	Height	Neck Circumference
00-1737-002-00	Small	3 in. (8cm)	16 1/2 in. (42cm)
00-1737-003-00	Medium	4 in. (10cm)	20 1/2 in. (52cm)
00-1737-004-00	Large	5 in. (13cm)	21 1/2 in. (55cm)
00-1737-001-00	Small, X-Wide	4 in. (10cm)	18 in. (46cm)
00-1737-005-00	Medium, X-Long	3 1/2 in. (9cm)	23 1/2 in. (60cm)

MEDIUM DENSITY			
Prod. No.	Size	Height	Neck Circumference
00-0737-002-00	Small	3 in. (8cm)	16 1/2 in. (42cm)
00-0737-003-00	Medium	4 in. (10cm)	20 1/2 in. (52cm)
00-0737-004-00	Large	5 in. (13cm)	21 1/2 in. (55cm)
00-0737-001-00	Small, X-Wide	4 in. (10cm)	18 in. (46cm)
00-0737-005-00	Universal	3 1/2 in. (9cm)	23 1/2 in. (60cm)
00-0737-006-00	Small, X-Long	3 in. (8cm)	21 in. (53cm)
00-0737-007-00	Low Contour	2 1/2 in. (6cm)	16 1/2 in. (42cm)



CERVICAL COLLARS

ZIMMER® STRAIGHT DELUXE FOAM COLLAR

- Maintains head and neck in neutral position
- Choice of 2 foam densities – firm or medium
- Designed to provide immobilization and gentle support of the cervical spine
- Washable
- Convenient contact closures for easy adjustment and application
- Indicated for cervical facet syndrome, muscle sprains and strains, degenerative disc disease with and without radiculopathy, discogenic pain



FIRM DENSITY			
Prod. No.	Size	Height	Neck Circumference
00-2740-002-00	Small	3 in. (8cm)	15 1/2 in. (39cm)
00-2740-003-00	Medium	4 in. (10cm)	20 1/2 in. (52cm)
00-2740-004-00	Large	5 in. (13cm)	21 1/2 in. (55cm)
00-2740-001-00	Small, X-Wide	4 in. (10cm)	18 in. (46cm)
00-2740-005-00	Medium, X-Long	3 1/2 in. (9cm)	23 1/2 in. (60cm)

MEDIUM DENSITY			
Prod. No.	Size	Height	Neck Circumference
00-2739-002-00	Small	3 in. (8cm)	15 1/2 in. (39cm)
00-2739-003-00	Medium	4 in. (10cm)	20 1/2 in. (52cm)
00-2739-004-00	Large	5 in. (13cm)	21 1/2 in. (55cm)
00-2739-001-00	Small, X-Wide	4 in. (10cm)	18 in. (46cm)
00-2739-005-00	Medium, X-Long	3 1/2 in. (9cm)	23 1/2 in. (60cm)

UNIVERSAL STRAIGHT FOAM COLLAR

- Provides mild to moderate support to relieve muscle strains and sprains
- Medium density foam
- Universal adult size for reduced inventories
- Collars packaged with extra cotton stockinette and extension piece



MEDIUM DENSITY			
Prod. No.	Size	Height	Neck Circumference
00-8705-045-00	Universal	3 1/2 in. (9cm)	21 1/2 in. (55cm)

PEDIATRIC STRAIGHT FOAM COLLAR

- Provides mild to moderate support to relieve muscle strains and sprains
- Medium density foam
- Universal pediatric size for reduced inventories



MEDIUM DENSITY			
Prod. No.	Size	Height	Neck Circumference
00-8601-380-00	Universal	2 in. (5cm)	12 in. (30cm)

UPPER EXTREMITY/SHOULDER/CLAVICLE SUPPORTS

ZIMMER® UNIVERSAL MULTI-SLING

- Envelope adjusts with contact closures for quick, easy universal sizing that reduces the need to stock multiple sizes in inventory
- Adjustable strap with contact closures for easy application and positioning
- Lightweight, breathable blue Tietex material aids air circulation to help keep the arm cool and comfortable
- Entire strap is Velfoam padded for patient comfort and compliance
- May be custom-imprinted with your logo as a special



Prod. No.	Size	Width	Envelope
00-1798-001-00	Universal	8 in. (20cm)	11 3/8 in. (29cm) Minimum and 19 in. (48cm) Maximum

AMERICAN WAY ARM SLING

- Envelope adjusts with contact closures for quick, easy universal sizing that reduces the need to stock multiple sizes in inventory
- Adjustable strap with contact closures for easy application and positioning
- Lightweight, breathable patriotic Tietex material aids air circulation to help keep the arm cool and comfortable
- Entire strap is Velfoam padded for improved patient comfort and compliance



Prod. No.	Size	Width	Envelope
00-1898-001-00	Universal	8 in. (20cm)	11 3/8 in. (29cm) Minimum and 19 in. (48cm) Maximum

ZIMMER® PADDED STRAP ARM SLING

- Deep envelope and adjustable strap
- Thumb loop inside envelope minimizes sling migration
- Contact closure for easy application
- Constructed of light-weight, breathable blue Tietex material which aids air circulation
- Entire strap is Velfoam padded for patient comfort and compliance



Prod. No.	Size	Envelope
00-1892-001-00	Small	6 in. x 13 in. (15cm x 33cm)
00-1892-002-00	Medium	7 in. x 16 in. (18cm x 41cm)
00-1892-003-00	Large	8 in. x 19 in. (20cm x 48cm)
00-1892-004-00	X-Large	9 in. x 22 in. (23cm x 56cm)



UPPER EXTREMITY/SHOULDER/CLAVICLE SUPPORTS

FASHION ARM SLING

- Thumb loop inside envelope minimizes sling migration
- Deep pocket for comfort and fit
- Lightweight, breathable blue Tietex material aids air circulation to keep arm cool and comfortable
- Also available in convenient 10-pack
- Contact closure
- Wide 1 1/2 in. (4cm) strap webbing helps distribute weight evenly for patient comfort and fit
- May be custom-imprinted with your logo as a special ✧



Prod. No.		Size	Envelope
Individual	Box of 10		
00-1796-011-00	00-1796-011-10	Small	6 in. x 13 in. (15cm x 33cm)
00-1796-012-00	00-1796-012-10	Medium	7 in. x 16 in. (18cm x 41cm)
00-1796-013-00	00-1796-013-10	Large	8 in. x 19 in. (20cm x 48cm)
00-1796-014-00	00-1796-014-10	X-Large	9 in. x 22 in. (23cm x 56cm)

PEDIATRIC FASHION ARM SLING

- Wide webbing shoulder strap for comfort and fit
- Easy slide buckle adjustment and contact closure allows wide range of lengths
- Available in popular child's print



Prod. No.	Color	Size	Envelope
00-8601-377-00	Navy	X-Small	4 1/2 in. x 9 1/2 in. (11cm x 24cm)
00-8601-377-50	Child's Print	X-Small	4 1/2 in. x 9 1/2 in. (11cm x 24cm)

ZIMMER® ARM SLING

- Single 1 in. (3cm) web strap crosses the back and fastens on envelope above the elbow for better weight distribution
- Lightweight, breathable Tietex material aids air circulation to help keep arm cool and comfortable
- Sold in boxes of 10
- May be custom-imprinted with your logo as a special ✧



Prod. No.	Size	Envelope
00-0764-001-10	Small	7 1/2 in. x 13 1/2 in. (19cm x 34cm)
00-0764-002-10	Medium	8 1/4 in. x 16 1/2 in. (21cm x 42cm)
00-0764-003-10	Large	9 1/2 in. x 20 in. (24cm x 51cm)
00-0764-004-10	X-Large	10 1/4 in. x 23 1/2 in. (26cm x 60cm)

UPPER EXTREMITY/SHOULDER/CLAVICLE SUPPORTS

ZIMMER® PEDIATRIC ARM SLING

- Single 1 in. (3cm) web strap crosses the back and fastens on envelope above the elbow for better weight distribution
- Natural color, lightweight 50/50 poly cotton material for easy care



Prod. No.	Size	Envelope
00-0764-010-00	XX-Small, Pediatric	5 1/2 in. x 8 in. (14cm x 20cm)
00-0764-011-00	X-Small, Pediatric	6 1/2 in. x 10 in. (17cm x 25cm)

STRAP ARM SLING WITH NECK PAD

- Includes perforated vinyl neck pad
- Fully adjustable strap arm sling
- Ideal for cast support



Prod. No.	Size	Strap Width
66-0576-001-00	Small	1 1/4 in. (3cm)
66-0576-002-00	Medium	1 1/2 in. (4cm)
66-0576-003-00	Large	1 1/2 in. (4cm)
66-0576-004-00 ◇	X-Large	2 in. (5cm)

ZIMFOAM® TUBE STRAP PADS

- Foam tube fits over the neck strap of an arm sling for added cushion and patient comfort
- Slotted tube makes it easy to apply and adjust to the appropriate position on the strap
- Packaged 4 per bag



Prod. No.	Size	Length
00-1086-022-00	Small	10 in. (25cm)
00-1086-023-00	Medium	12 in. (25cm)
00-1086-024-00	Large	15 in. (38cm)

ZIMMER® PADDED SHOULDER IMMOBILIZER WITH PADDED BODY STRAP

- Resists abduction and external rotation
- Lightweight, breathable Tietex material aids air circulation to help keep arm cool and comfortable
- Entire neck and body straps are Velfoam padded for patient comfort



Prod. No.	Size	Envelope
00-1893-001-00	Small	6 in. x 13 in. (15cm x 33cm)
00-1893-002-00	Medium	7 in. x 16 in. (18cm x 41cm)
00-1893-003-00	Large	8 in. x 19 in. (20cm x 48cm)
00-1893-004-00	X-Large	9 in. x 22 in. (23cm x 56cm)



UPPER EXTREMITY/SHOULDER/CLAVICLE SUPPORTS

SHOULDER IMMOBILIZER WITH BODY STRAP

- Resists abduction and external rotation
- Lightweight, breathable Tietex material aids air circulation to help keep arm cool and comfortable



Prod. No.	Size	Envelope
00-0763-001-00	Small	7 1/2 in. x 12 1/2 in. (19cm x 32cm)
00-0763-002-00	Medium	8 1/4 in. x 15 1/2 in. (21cm x 39cm)
00-0763-003-00	Large	9 1/2 in. x 19 in. (24cm x 48cm)
00-0763-004-00	X-Large	10 in. x 22 1/2 in. (25cm x 57cm)
00-0763-011-00	X-Small, Pediatric	6 1/2 in. x 10 in. (17cm x 25cm)

DELUXE SLING AND SWATHE SHOULDER IMMOBILIZER

- Holds arm securely against body for shoulder immobilization
- Perforated foam allows air circulation for patient comfort
- Contact closures for easy adjustment
- Optional use as cradle sling



Prod. No.	Size	Envelope	Swathe
00-3898-005-00	Small	9 in. x 11 in. (23cm x 28cm)	4 in. x 32 in. (10cm x 81cm)
00-3898-006-00	Medium	9 1/2 in. x 11 3/4 in. (24cm x 30cm)	5 in. x 38 in. (13cm x 97cm)
00-3898-007-00	Large	10 in. x 13 in. (25cm x 33cm)	5 1/2 in. x 46 in. (14cm x 117cm)
00-3898-008-00	X-Large	10 in. x 13 in. (25cm x 33cm)	5 1/2 in. x 57 in. (14cm x 145cm)

PEDIATRIC DELUXE SLING AND SWATHE

- Holds arm securely against body for shoulder immobilization
- Perforated foam allows air circulation for patient comfort
- Contact closures for easy adjustment
- Optional use as cradle sling
- Washable



Prod. No.	Size	Envelope	Swathe
00-8601-796-00	X-Small	8 in. x 10 in. (20cm x 25cm)	4 in. x 32 in. (10cm x 81cm)

UNIVERSAL SLING AND SWATHE

- Perforated foam increases air circulation for patient comfort
- May be trimmed to fit patient
- Universal sizing to minimize inventory



Prod. No.	Size	Envelope	Swathe
00-3899-001-00	Universal	10 in. x 13 in. (25cm x 33cm)	5 in. x 57 in. (13cm x 145cm)

Prod. No.	Description
00-8614-001-71	Swath Only

UPPER EXTREMITY/SHOULDER/CLAVICLE SUPPORTS

ELASTIC SHOULDER IMMOBILIZER

- Effective support without neck strap for shoulder dislocations
- Chestband contoured for women



MEN'S		
Prod. No.	Size	Chest Range
00-0783-001-00	Small	24 in. - 30 in. (61cm - 76cm)
00-0783-002-00	Medium	30 in. - 36 in. (76cm - 91cm)
00-0783-003-00	Large	36 in. - 42 in. (91cm - 107cm)
00-0783-004-00	X-Large	42 in. - 48 in. (107cm - 122cm)
WOMEN'S		
Prod. No.	Size	Chest Range
00-0784-001-00	Small	24 in. - 30 in. (61cm - 76cm)
00-0784-002-00	Medium	30 in. - 36 in. (76cm - 91cm)
00-0784-003-00	Large	36 in. - 42 in. (91cm - 107cm)
00-0784-004-00	X-Large	42 in. - 48 in. (107cm - 122cm)
00-0784-005-00 ◇	XX-Large	48 in. - 54 in. (122cm - 137cm)

UNIVERSAL FOAM SHOULDER IMMOBILIZER

- Effective support without neck strap for shoulder dislocations
- Wrist envelope securely attached to chest swathe
- Foam padding and removable shoulder straps enhance comfort



Prod. No.	Size	Chest Range
00-1785-001-00	Universal	24 in. - 46 in. (61cm - 117cm)

ZIMFOAM® SHOULDER ABDUCTION POSITIONER

- Adjust for 40 or 90 degrees
- Wide comfortable straps
- Can be trimmed if necessary



Prod. No.	Size
00-2675-000-00	Universal

GARDNER ARM ELEVATOR

- Helps control edema following hand trauma or surgery
- Universal size
- Contact closure facilitates application and arm inspection
- Standard version features lightweight, breathable Tietex material which aids in air circulation
- Deluxe version features canvas construction and includes attached ice pockets



Prod. No.	Size	Material	Envelope
00-8600-194-00	Ped	Tietex	6 in. x 10 in. (15cm x 25cm)
00-8600-195-00	Youth	Tietex	8 in. x 13 5/8 in. (20cm x 35cm)
00-1790-001-00 ◇	Adult, Standard	Tietex	8 1/4 in. x 14 1/2 in. (21cm x 37cm)
00-1790-002-00 ◇	Adult, Deluxe	Canvas	8 1/4 in. x 14 1/2 in. (21cm x 37cm)



UPPER EXTREMITY/SHOULDER/CLAVICLE SUPPORTS

UNIVERSAL FOAM ARM ELEVATOR/SLING

- Helps control edema following hand trauma or surgery
- Universal size
- Contact closure strap facilitates conversion from arm elevator to sling
- Large envelope can accommodate icebag
- Perforated foam increases air circulation



Prod. No.	Envelope
00-1978-001-00	10 1/4 in. x 19 1/2 in. (26cm x 50cm)

MCLEOD® PADDED CLAVICLE SPLINT

- Felt padding provides proper immobilization and enhances patient comfort
- Choice of buckle or contact closures
- Washable
- Indicated for most clavicle fractures without articular involvement and where ligaments are intact



Prod. No.				
Buckle	Contact Closures	Size	Strap Width	Strap Length
00-8600-282-00		XX-Small	1 3/4 in. (4cm)	15 1/2 in. (39cm)
00-0735-001-00	00-0735-011-00	X-Small (Infant)	1 5/8 in. (4cm)	20 1/2 in. (52cm)
00-0735-002-00	00-0735-012-00	Small (Child)	2 in. (5cm)	23 3/4 in. (60cm)
00-0735-003-00	00-0735-013-00	Medium (Youth)	2 1/4 in. (6cm)	27 1/4 in. (69cm)
00-0735-004-00	00-0735-014-00	Large (Adult)	2 3/4 in. (7cm)	34 1/2 in. (88cm)
00-0735-005-00	00-0735-015-00	X-Large	3 in. (8cm)	41 in. (104cm)
00-0735-006-00	◇	XX-Large	3 1/4 in. (8cm)	43 in. (109cm)

ZIMMER® PREMIUM GEL TENNIS ELBOW SUPPORT

- Gel pad for hot/cold therapy and added compression on the extensor without cutting off circulation
- Sporty black color
- Universal 2 in. (5cm) wide support with contact closure for easy application and perfect adjustment
- Washable for better patient hygiene



Prod. No.	Size
00-1445-001-00	Universal

UPPER EXTREMITY/SHOULDER/CLAVICLE SUPPORTS

ZIMMER® FOREARM STRAP (TENNIS ELBOW SUPPORT)

- Curvilinear design for accurate fit
- Wide band disseminates force, relieves extensor tension and pain
- Contact closure for proper fit and easy adjustment
- 2 3/4 in. (7cm) wide
- Foam with tricot lining
- Hand washable for better patient hygiene
- White color
- Sold individually or in boxes of 6



Prod. No.		Size	Forearm Circumference
Individual	Box of 6		
00-1295-008-00	00-1295-002-00	Small	8 1/2 in. - 9 1/2 in. (22cm - 24cm)
00-1295-009-00	00-1295-003-00	Medium	9 1/2 in. - 10 1/2 in. (24cm - 27cm)
00-1295-010-00	00-1295-004-00	Large	10 1/2 in. - 11 1/2 in. (27cm - 29cm)
00-1295-011-00	00-1295-005-00	X-Large	11 1/2 in. - 12 1/2 in. (29cm - 32cm)

FROIMSON TENNIS ELBOW SUPPORT

- 2 1/2 in. (6cm) wide band
- Natural color
- Hand washable for better patient hygiene
- Universal size minimizes inventory
- Contact closure for proper fit and easy adjustment
- Sold in boxes of 6



Prod. No.	Size	Forearm Circumference
00-0754-006-00	Universal	9 1/2 in. - 14 1/2 in. (24cm - 37cm)

NEOPRENE ELBOW SUPPORT

- 1/8 in. (3mm) thick neoprene provides compression and warmth for mild to moderate elbow strains and sprains
- Circumferential 2 in. (5cm) strap
- Contact closure provides exact fit and adjustment
- 11 in. (28cm) long pull-on sleeve
- Quality construction for superior durability
- Anatomically contoured fit
- Beige color



Prod. No.	Size	Forearm Circumference
00-3448-003-00	Small	9 in. - 10 in. (23cm - 25cm)
00-3448-004-00	Medium	10 in. - 11 in. (25cm - 28cm)
00-3448-005-00	Large	11 in. - 12 in. (28cm - 30cm)
00-3448-006-00	X-Large	12 in. - 14 in. (30cm - 36cm)



WRIST SUPPORTS

ZIMMER® PREMIUM WRIST SPLINT


- Constructed of a lightweight and durable fabric with a soft lycra-lined interior that is extremely breathable
- 3-dimensional molding technology varies padding thickness for superior comfort and fit
- Adjustable design for customized fit
- Contoured along distal palmar crease to allow for 90 degrees MP flexion
- Precisely sculpted narrow web space helps prevent pinching
- V-notch at 5th MP to allow for palmar spread
- Conformable radial and adjustable palmar stays
- Contact closure straps for easy application and removal



Prod. No.			
Left	Right	Size	Wrist Circumference
00-1476-001-00	00-1476-002-00	X-Small	5 1/4 in. - 6 1/4 in. (13cm - 16cm)
00-1476-003-00	00-1476-004-00	Small	6 1/4 in. - 7 1/4 in. (16cm - 18cm)
00-1476-005-00	00-1476-006-00	Medium	7 1/4 in. - 8 in. (18cm - 20cm)
00-1476-007-00	00-1476-008-00	Large	8 in. - 9 in. (20cm - 23cm)
00-1476-009-00	00-1476-010-00	X-Large	9 in. - 10 in. (23cm - 25cm)

ZIMMER® PREMIUM THUMB SPICA

- Constructed of a lightweight and durable fabric with a soft lycra-lined interior that is extremely breathable
- 3-dimensional molding technology varies padding thickness for exceptional comfort and fit
- Adjustable design for customized fit
- Indicated for Gamekeeper's thumb, DeQuervain's tendonitis, selected scaphoid injuries, post cast healing and soft tissue injuries
- Designed for effective long-term immobilization and protection
- Conformable radial and adjustable palmar stays
- Contact closure straps for easy application and removal



Prod. No.			
Left	Right	Size	Wrist Circumference
00-1477-001-00	00-1477-002-00	X-Small	5 1/4 in. - 6 1/4 in. (13cm - 16cm)
00-1477-003-00	00-1477-004-00	Small	6 1/4 in. - 7 1/4 in. (16cm - 18cm)
00-1477-005-00	00-1477-006-00	Medium	7 1/4 in. - 8 in. (18cm - 20cm)
00-1477-007-00	00-1477-008-00	Large	8 in. - 9 in. (20cm - 23cm)
00-1477-009-00	00-1477-010-00	X-Large	9 in. - 10 in. (23cm - 25cm)

WRIST SUPPORTS

ZIMMER® PREMIUM EXO-WRIST BRACE

- Indicated for the treatment of strains, sprains and repetitive stress injuries, including Carpal Tunnel Syndrome
- Lightweight low profile design helps reduce bulk and promotes patient compliance
- Pre-shaped to maintain wrist in a neutral position, important for treatment
- Exo-skeleton is contoured to ensure anatomically correct fit and comfortable limited motion and flexibility resisting extremes in range of motion
- Exo-skeleton provides radial and ulnar styloid pressure relief while providing more room at the base of the thumb and palmar crease
- 3-dimensional molding technology varies padding thickness for superior comfort and fit
- Web space is thin and soft to help reduce pinching and irritation
- Attractive appearance for better patient compliance
- Brace and removable padding are easily cleaned and maintained by using mild soap and warm water



Prod. No.		Size	Palm Circumference	Wrist Circumference
Left	Right			
00-1478-001-00	00-1478-002-00	Small	6 3/4 in. - 8 in. (17cm - 20cm)	5 1/4 in. - 6 1/2 in. (13cm - 17cm)
00-1478-003-00	00-1478-004-00	Medium	8 in. - 9 1/4 in. (20cm - 23cm)	6 1/2 in. - 7 3/4 in. (17cm - 20cm)
00-1478-005-00	00-1478-006-00	Large	9 1/4 in. - 10 1/2 in. (23cm - 27cm)	7 3/4 in. - 10 in. (20cm - 25cm)

ADJUSTABLE WRIST BRACE

- Immobilizes wrist and carpal joints while providing uniform support
- Malleable, removable metal stay reverses to "cocked-down" position
- Contoured to provide uniform support allowing full finger motion
- Elastic tension varies to provide maximum support to wrist without pinching hand
- Contact closures for easy application and adjustment
- Hand washable with mild soap and water



Prod. No.		Size	Wrist Circumference
Left	Right		
00-1799-001-00	00-1799-002-00	X-Small, Pediatric	5 in. - 6 in. (13cm - 15cm)
00-1799-003-00	00-1799-004-00	Small	5 1/2 in. - 6 1/2 in. (14cm - 17cm)
00-1799-005-00	00-1799-006-00	Medium	6 3/4 in. - 7 3/4 in. (17cm - 20cm)
00-1799-007-00	00-1799-008-00	Large	7 1/2 in. - 8 1/2 in. (19cm - 22cm)
00-1799-009-00	00-1799-010-00	X-Large	8 in. - 9 in. (20cm - 23cm)



WRIST SUPPORTS

UNIVERSAL COCKUP WRIST SPLINT

- 3 forearm straps for greater support and security
- Perforations provide for greater air circulation for patient comfort and compliance
- Easily applied with 1 hand
- Palmar and dorsal aluminum stays for maximum support



ADULT

Prod. No.		Size
Left	Right	
00-1769-003-00	00-1769-002-00	Universal

PEDIATRIC

Prod. No.		Size
Left	Right	
00-8601-887-03	00-8601-887-02	Universal

ZIMMER® COCKUP WRIST SPLINT

- 1 1/2 in. (4cm) aluminum stays hold wrist in dorsiflexion
- Soft felt material for comfort
- Aluminum stay for maximum support



Prod. No.

Left	Right	Size	Wrist Circumference
00-0799-001-00	00-0799-002-00	Small	5 in. - 7 1/2 in. (13cm - 19cm)
00-0799-003-00	00-0799-004-00	Medium	6 1/2 in. - 8 1/2 in. (17cm - 22cm)
00-0799-005-00	00-0799-006-00	Large	7 1/2 in. - 9 1/2 in. (19cm - 24cm)
00-0799-007-00	00-0799-008-00	X-Large	9 in. - 11 in. (23cm - 28cm)

DELUXE COCKUP WRIST SPLINT

- 1 1/2 in. (4cm) contoured aluminum stays hold wrist in dorsiflexion, yet permit natural extension of fingers
- Perforated vinyl allows increased air circulation for patient comfort and compliance
- 3 contact closures for firm immobilization



ADULT

Prod. No.		Size	Wrist Circumference
Left	Right		
00-1759-001-00	00-1759-002-00	Small	5 in. - 7 1/2 in. (13cm - 19cm)
00-1759-003-00	00-1759-004-00	Medium	6 1/2 in. - 8 1/2 in. (17cm - 22cm)
00-1759-005-00	00-1759-006-00	Large	7 1/2 in. - 9 1/2 in. (19cm - 24cm)
00-1759-007-00	00-1759-008-00	X-Large	9 in. - 11 in. (23cm - 28cm)

PEDIATRIC

Prod. No.		Size
Left	Right	
00-8600-529-00	00-8600-530-00	X-Small

WRISTLET WITH STRAP

- 1/8 in. (3mm) thick neoprene
- 4 contact closures
- 2 in. (5cm) contact strap for additional support
- Beige color



Prod. No.

Left	Right	Size	Wrist Circumference
00-3447-403-00	00-3447-303-00	Small	5 in. - 6 in. (13cm - 15cm)
00-3447-404-00	00-3447-304-00	Medium	6 in. - 7 in. (15cm - 18cm)
00-3447-405-00	00-3447-305-00	Large	7 in. - 8 in. (18cm - 20cm)
00-3447-406-00	00-3447-306-00	X-Large	8 in. - 9 in. (20cm - 23cm)

ABDOMINAL/BACK/RIB SUPPORTS

ZIMMER® COMFORT PANEL RIB BELT

- Universal size, fits 28 in. to 57 in. chests (71cm - 145cm)
- 2-panel special elastic construction provides even support and compression
- Women's rib belts are contoured
- 6 in. (15cm) wide with a comfort panel lining for a comfortable, secure fit



MEN'S			
Prod. No.	Size	Width	Chest Range
00-0584-011-00	Universal	6 in. (15cm)	28 in. - 57 in. (71cm - 145cm)
WOMEN'S			
Prod. No.	Size	Width	Chest Range
00-0584-012-00	Universal	6 in. (15cm)	28 in. - 57 in. (71cm - 145cm)

UNIVERSAL RIB BELT

- 6 vertical pile strips for maximum circumference adjustment
- Available in 6 in. and 8 in. widths (15cm and 20cm)
- Universal size
- Packaged in boxes of 6



MEN'S			
Prod. No.	Size	Width	Chest Range
00-2721-006-00	Universal	6 in. (15cm)	24 in. - 48 in. (61cm - 122cm)
00-2721-008-00	Universal	8 in. (20cm)	24 in. - 48 in. (61cm - 122cm)
WOMEN'S			
Prod. No.	Size	Width	Chest Range
00-2722-006-00	Universal	6 in. (15cm)	24 in. - 48 in. (61cm - 122cm)
00-2722-008-00	Universal	8 in. (20cm)	24 in. - 48 in. (61cm - 122cm)

ZIMMER® COMFORT PANEL ABDOMINAL SUPPORT

- Special elastic panel construction provides even support
- Available in 9 in. (23cm) width with 3-panel construction or 12 in. (30cm) width with 4-panel construction
- Contact closure adjustment possible at any point on 12 in. (30cm) long comfort panel



3-PANEL CONSTRUCTION			
Prod. No.	Size	Width	Length
00-0584-031-00	Medium	9 in. (23cm)	30 in. - 45 in. (76cm - 114cm)
00-0584-032-00	Large	9 in. (23cm)	46 in. - 62 in. (117cm - 157cm)
4-PANEL CONSTRUCTION			
Prod. No.	Size	Width	Length
00-0584-041-00	Medium	12 in. (30cm)	30 in. - 45 in. (76cm - 114cm)
00-0584-042-00	Large	12 in. (30cm)	46 in. - 62 in. (117cm - 157cm)
00-0584-043-00 ◊	X-Large	12 in. (30cm)	62 in. - 77 in. (157cm - 196cm)
00-0584-044-00 ◊	XX-Large	12 in. (30cm)	77 in. - 94 in. (196cm - 239cm)

Measurements taken at waist for men and at hips for women

UNIVERSAL ABDOMINAL SUPPORT

- Elastic construction assures firm fit
- Velfoam tabs for maximum sizing adjustments
- Universal size minimizes inventory
- Available in 10 in. (25cm) or 12 in. (30cm) widths
- Packaged in boxes of 6



Prod. No.	Size	Width	Length
00-2732-010-00	Universal	10 in. (25cm)	26 in. - 50 in. (66cm - 127cm)
00-2732-012-00	Universal	12 in. (30cm)	26 in. - 50 in. (66cm - 127cm)

Measurements taken at waist for men and at hips for women



ABDOMINAL/BACK/RIB SUPPORTS

ELASTIC ABDOMINAL SACRAL SUPPORT

- Available in 10 in. (25cm) width only
- Stays prevent edge rollover
- Contact closures for easy application
- Elastic construction enhances firm, comfortable fit



Prod. No.	Width	Length
00-0734-000-28 ◇	10 in. (25cm)	28 in. - 32 in. (71cm - 81cm)
00-0734-000-32 ◇	10 in. (25cm)	32 in. - 36 in. (81cm - 91cm)
00-0734-000-36 ◇	10 in. (25cm)	36 in. - 40 in. (91cm - 102cm)
00-0734-000-40 ◇	10 in. (25cm)	40 in. - 45 in. (102cm - 114cm)
00-0734-000-45 ◇	10 in. (25cm)	45 in. - 50 in. (114cm - 127cm)

ZIMMER® PREMIUM INFLATABLE BACK SUPPORT WITH GEL PAD

- Adjustable air bladder provides a custom fit and offers focused compression in the lumbar area
- Self-contained pneumatic system designed to be low profile and easy to use
- Gel pad for hot and cold therapy is included
- Mesh panel and ventilation holes increase breathability and enhance comfort
- Overlapping double pull strap design provides additional support to the lumbar and abdominal areas
- Rubber "gripper strips" sewn inside each support help prevent migration



Prod. No.	Description	Waist
00-1479-001-00	Small	25 in. - 30 in. (64cm - 76cm)
00-1479-002-00	Medium	30 in. - 34 in. (76cm - 86cm)
00-1479-003-00	Large	34 in. - 38 in. (86cm - 97cm)
00-1479-004-00	X-Large	38 in. - 42 in. (97cm - 107cm)
00-1479-005-00	XX-Large	42 in. - 56 in. (107cm - 142cm)

UNIVERSAL SACRO-LUMBAR SUPPORT

- Provides lower back support for patients with up to 48 in. (122cm) hip measurement
- Contoured movable stays for custom fit and support
- 11 in. (28cm wide)
- Converts to pelvic belt by attaching Traction Strap Kit



Prod. No.	Description	Size	Width	Length
00-1730-001-00	Sacro-Lumbar Support	Universal	11 in. (28cm)	Up to 48 in. (122cm)
00-1570-011-00	Traction Strap Kit (Box 2)			

ABDOMINAL/BACK/RIB SUPPORTS

SPS III

- 4 rigid contoured steel stays for maximum support
- 13 in. (33 cm) width
- 3 adjustment straps with contact closures for secure fit



Prod. No.	Size	Measurements
00-1735-021-00	X-Small	28 in. - 32 in. (71cm - 81cm)
00-1735-022-00	Small	32 in. - 36 in. (81cm - 91cm)
00-1735-023-00	Medium	36 in. - 40 in. (91cm - 102cm)
00-1735-024-00	Large	40 in. - 44 in. (102cm - 112cm)
00-1735-025-00	X-Large	44 in. - 48 in. (112cm - 122cm)
00-1735-026-00 ◊	XX-Large	48 in. - 52 in. (122cm - 132cm)
00-1735-027-00 ◊	XXX-Large	52 in. - 56 in. (132cm - 142cm)
00-1735-028-00 ◊	XXXX-Large	56 in. - 60 in. (142cm - 152cm)

Measurements taken at waist for men and at hips for women

ELASTIC BACK SUPPORT

- Side panel elastic adjustment straps provide uniform compression
- 13 in. (33cm) width
- 2 flexible steel stays
- Convenient contact closures for easy application



Prod. No.	Size	Measurements
00-2735-021-00 ◊	X-Small	28 in. - 32 in. (71cm - 81cm)
00-2735-022-00	Small	32 in. - 36 in. (81cm - 91cm)
00-2735-023-00	Medium	36 in. - 40 in. (91cm - 102cm)
00-2735-024-00	Large	40 in. - 44 in. (102cm - 112cm)
00-2735-025-00	X-Large	44 in. - 48 in. (112cm - 122cm)

Measurements taken at waist for men and at hips for women

UNIVERSAL LUMBAR SUPPORT

- Indicated for facet syndrome, muscle spasm, stiffness, low back pain, or post-operative support and compression
- 1/8 in. (3mm) thick neoprene
- 1/2 in. (13mm) sewn-in lumbar pad for additional support
- Beige color



Prod. No.	Size	Measurements
00-3449-001-00	Small/Medium	20 in. - 40 in. (51cm - 102cm)
00-3449-002-00	Large/X-Large	22 in. - 48 in. (56cm - 122cm)

Measurements taken at waist for men and at hips for women



KNEE IMMOBILIZERS

PERFORATED FOAM-LINED KNEE IMMOBILIZER

- Breathable perforated foam allows increased air circulation for added patient comfort
- Lightweight construction
- Radiolucent aluminum stays
- Available with contoured or straight posterior stays
- Leg circumference sizes range from Small to XX-Large



Size	Thigh Circumference
Small	18 in. - 20 in. (46cm - 51cm)
Medium	20 in. - 22 in. (51cm - 56cm)
Large	22 in. - 24 in. (56cm - 61cm)
X-Large	24 in. - 26 in. (61cm - 66cm)
XX-Large	26 in. - 28 in. (66cm - 71cm)

CONTOURED STAYS

Size	12 in. (30cm) Length	16 in. (41cm) Length	19 in. (48cm) Length	22 in. (56cm) Length	24 in. (61cm) Length
Small	00-1876-503-12	00-1876-503-16	00-1876-503-19	00-1876-503-22	00-1876-503-24
Medium	00-1876-504-12	00-1876-504-16	00-1876-504-19	00-1876-504-22	00-1876-504-24
Large	00-1876-505-12	00-1876-505-16	00-1876-505-19	00-1876-505-22	00-1876-505-24
X-Large	00-1876-506-12	00-1876-506-16	00-1876-506-19	00-1876-506-22	00-1876-506-24
XX-Large			00-1876-507-19	00-1876-507-22	

STRAIGHT STAYS

Size	12 in. (30cm) Length	16 in. (41cm) Length	19 in. (48cm) Length	22 in. (56cm) Length	24 in. (61cm) Length
Small	00-1876-603-12		00-1876-603-19	00-1876-603-22	00-1876-603-24
Medium	00-1876-604-12	00-1876-604-16	00-1876-604-19	00-1876-604-22	00-1876-604-24
Large	00-1876-605-12	00-1876-605-16	00-1876-605-19	00-1876-605-22	00-1876-605-24
X-Large	00-1876-606-12		00-1876-606-19	00-1876-606-22	00-1876-606-24
XX-Large			00-1876-607-19		

TRICOT-LINED KNEE IMMOBILIZER

- Breathable foam with inner lining of soft nylon tricot
- Foam interior medial/lateral strips help keep immobilizer correctly positioned and help to eliminate distal migration
- Radiolucent contoured aluminum posterior stays keep leg in comfortable neutral position
- Leg circumference sizes range from Small to XX-Large



Size	Thigh Circumference
Small	18 in. - 20 in. (46cm - 51cm)
Medium	20 in. - 22 in. (51cm - 56cm)
Large	22 in. - 24 in. (56cm - 61cm)
X-Large	24 in. - 26 in. (61cm - 66cm)
XX-Large	26 in. - 28 in. (66cm - 71cm)

CONTOURED STAYS

Size	12 in. (30cm) Length	16 in. (41cm) Length	19 in. (48cm) Length	22 in. (56cm) Length	24 in. (61cm) Length
Small	00-1877-503-12	00-1877-503-16	00-1877-503-19	00-1877-503-22	00-1877-503-24
Medium	00-1877-504-12	00-1877-504-16	00-1877-504-19	00-1877-504-22	00-1877-504-24
Large	00-1877-505-12	00-1877-505-16	00-1877-505-19	00-1877-505-22	00-1877-505-24
X-Large	00-1877-506-12	00-1877-506-16	00-1877-506-19	00-1877-506-22	00-1877-506-24
XX-Large			00-1877-507-19	00-1877-507-22	

KNEE IMMOBILIZERS

UNIVERSAL KNEE IMMOBILIZER WITH MOVABLE MEDIAL/LATERAL STAYS

- Lightweight, breathable foam/nylon material
- Three medial/lateral stay pocket options
- Radiolucent contoured aluminum stays
- Contoured posterior stays keep leg in comfortable neutral position



Prod. No.	Length	Max. Thigh Circumference
00-1880-500-12	12 in. (30cm)	28 in. (71cm)
00-1880-500-16	16 in. (41cm)	28 in. (71cm)
00-1880-500-19	19 in. (48cm)	29 in. (74cm)
00-1880-500-22	22 in. (56cm)	29 in. (74cm)
00-1880-500-24	24 in. (61cm)	30 in. (76cm)
00-1880-500-26	26 in. (66cm)	30 in. (76cm)
00-1880-500-28	28 in. (71cm)	31 in. (79cm)

UNIVERSAL KNEE IMMOBILIZER WITH TRIMMABLE FOAM

- Breathable trimmable foam for patient comfort and custom fit
- Lightweight construction
- Radiolucent aluminum stays
- Contoured and straight posterior stays available



CONTOURED STAYS		
Prod. No.	Length	Max. Thigh Circumference
00-1882-500-12	12 in. (30cm)	29 in. (74cm)
00-1882-500-16	16 in. (41cm)	29 in. (74cm)
00-1882-500-19	19 in. (48cm)	30 in. (76cm)
00-1882-500-22	22 in. (56cm)	30 in. (76cm)
00-1882-500-24	24 in. (61cm)	30 in. (76cm)
00-1882-500-26	26 in. (66cm)	31 in. (79cm)
00-1882-500-28	28 in. (71cm)	32 in. (81cm)

STRAIGHT STAYS		
Prod. No.	Length	Max. Thigh Circumference
00-1882-600-12	12 in. (30cm)	29 in. (74cm)
00-1882-600-16	16 in. (41cm)	29 in. (74cm)
00-1882-600-19	19 in. (48cm)	30 in. (76cm)
00-1882-600-22	22 in. (56cm)	30 in. (76cm)
00-1882-600-24	24 in. (61cm)	30 in. (76cm)

UNIVERSAL PERFORATED TRI-PANEL KNEE IMMOBILIZER

- Adjustable perforated lightweight panels are designed to accommodate cone-shaped legs and peel back to allow quick and easy access to wound site
- Breathable perforated foam and perforated vinyl allows increased air circulation for added patient comfort
- Contoured posterior stays keep leg in comfortable neutral position
- Radiolucent aluminum stays
- Contact closures for easy application, adjustability and optimum fit



Prod. No.	Length	Max. Thigh Circumference
00-1883-500-12	12 in. (30cm)	27 in. (69cm)
00-1883-500-16	16 in. (41cm)	27 in. (69cm)
00-1883-500-19	19 in. (48cm)	27 in. (69cm)
00-1883-500-22	22 in. (56cm)	27 in. (69cm)
00-1883-500-24	24 in. (61cm)	28 in. (71cm)
00-1883-500-26	26 in. (66cm)	28 in. (71cm)
00-1883-500-28	28 in. (71cm)	28 in. (71cm)



KNEE IMMOBILIZERS

UNIVERSAL KNEE IMMOBILIZER WITH ELASTIC PATELLA STRAP

- Breathable, perforated foam/nylon material
- Elastic patella strap for secure positioning
- Radiolucent contoured aluminum posterior stays



Prod. No.	Length	Max. Thigh Circumference
00-1885-500-12	12 in. (30cm)	28 in. (71cm)
00-1885-500-16	16 in. (41cm)	28 in. (71cm)
00-1885-500-19	19 in. (48cm)	29 in. (74cm)
00-1885-500-22	22 in. (56cm)	29 in. (74cm)
00-1885-500-24	24 in. (61cm)	30 in. (76cm)
00-1885-500-26	26 in. (66cm)	30 in. (76cm)
00-1885-500-28	28 in. (71cm)	31 in. (79cm)

UNIVERSAL KNEE IMMOBILIZER WITH ELASTIC STRAPS

- 3 in. (8cm) wide elastic straps that maintain compression and support
- Perforated foam helps increase air circulation for patient comfort and compliance
- Absence of "D" rings make application quick and easy



Prod. No.	Length	Max. Thigh Circumference
00-1888-500-12	12 in. (30cm)	29 in. (74cm)
00-1888-500-16	16 in. (41cm)	29 in. (74cm)
00-1888-500-19	19 in. (48cm)	30 in. (76cm)
00-1888-500-22	22 in. (56cm)	30 in. (76cm)
00-1888-500-24	24 in. (61cm)	30 in. (76cm)
00-1888-500-26	26 in. (66cm)	31 in. (79cm)
00-1888-500-28	28 in. (71cm)	32 in. (81cm)

POPLITEAL FOAM INSERT

- Inserts in the popliteal area of the knee immobilizer to provide a better fit
- Helps prevent distal migration of the knee immobilizer
- Contoured shape for better fit
- Packaged 10 per box



Prod. No.
00-1889-000-10

KNEE BRACES AND SUPPORTS

DROP-LOK® KNEE BRACE

- Removable hinge system offers maximum stability in full extension when locked, lateral stability with limited flexion when unlocked
- Foam cartilage pad for added patellar support



Prod. No.	Size	Knee Circumference
00-1746-001-00	X-Small	9 1/2 in. - 11 1/2 in. (24cm - 29cm)
00-1746-002-00	Small	11 1/2 in. - 13 in. (29cm - 33cm)
00-1746-003-00	Medium	13 in. - 14 3/4 in. (33cm - 37cm)
00-1746-004-00	Large	14 3/4 in. - 16 1/2 in. (37cm - 42cm)
00-1746-005-00	X-Large	16 1/2 in. - 18 1/2 in. (42cm - 47cm)

Measured around the center of the knee cap with the leg fully extended

ELASTIC KNEE SUPPORT

- Elastic support with medial and lateral support straps
- Semi-circular patellar cut-out
- Contact closures for easy application



Prod. No.	Size	Knee Circumference
00-0774-001-00	Small	12 in. - 14 in. (30cm - 36cm)
00-0774-002-00	Medium	14 in. - 16 in. (36cm - 41cm)
00-0774-003-00	Large	16 in. - 18 in. (41cm - 46cm)
00-0774-004-00	X-Large	18 in. - 20 in. (46cm - 51cm)

Measured around the center of the knee cap with the leg fully extended

HINGED KNEE SUPPORT-HEAVY DUTY

- Heavy-duty steel hinges for added lateral stability
- Allows knee to tolerate more strenuous activity during recuperation
- Lace and strap closures



Prod. No.	Size	Knee Circumference
00-1773-002-00	Small	13 in. (33cm)
00-1773-003-00	Medium	14 in. (36cm)
00-1773-004-00	Large	15 1/2 in. (39cm)
00-1773-005-00	X-Large	16 1/2 in. (42cm)
00-1773-006-00 ⇄	XX-Large	17 1/2 in. (44cm)

Measured around the center of the knee cap with the leg fully extended

CARTILAGE KNEE BRACE

- Helical spring steel stays permit normal flexion with rotational stability
- Foam cartilage pad for added patellar support
- Indicated for minor injuries involving the menisci or ligaments



Prod. No.	Size	Knee Circumference
00-1747-001-00	X-Small, Pediatric	9 1/2 in. - 11 in. (24cm - 29cm)
00-1747-002-00	Small	11 1/2 in. - 13 in. (29cm - 33cm)
00-1747-003-00	Medium	13 in. - 14 3/4 in. (33cm - 37cm)
00-1747-004-00	Large	14 3/4 in. - 16 1/2 in. (37cm - 42cm)
00-1747-005-00	X-Large	16 1/2 in. - 18 1/2 in. (42cm - 47cm)

Measured around the center of the knee cap with the leg fully extended



KNEE BRACES AND SUPPORTS

POST-OP FLEXION/EXTENSION KNEE BRACE

- Post-operative brace effective for rehabilitation, joint stiffness, and contractures, instability, strains, sprains and ligament repairs
- Range of motion hinge with flexion/extension stop settings and static lock pin for immobilization or static progress
- Range of motion from 0 degrees to 110 degrees in 10 degree increments
- Universal size
- Telescopic length adjustment
- Lightweight design for patient compliance with plastic panels
- Soft, open-cell padding for patient comfort
- Strong, flexible paddles provide solid support



Prod. No.	Length	Maximum Thigh Circumference	Maximum Calf Circumference
00-1558-001-00	24 1/2 in. - 27 1/2 in. (62cm - 70cm)	27 in. (69cm)	20 in. (51cm)
00-1558-002-00	19 in. - 21 1/2 in. (48cm - 55cm)	27 in. (69cm)	20 in. (51cm)

ZIMMER® PROPOINT* KNEE BRACE

- Adjustable flexion and extension ROM control
- Internal contoured cuffs create added stability
- Accordion molding in the popliteal region eliminates bunching and provides exceptional patient comfort
- Fully integrated 1-piece design is easy to apply
- Indicated for medial or lateral collateral ligament instabilities, anterior or posterior cruciate ligament sprains, and meniscal injuries
- 3-dimensional molding technology varies padding thickness for superior comfort and fit
- Available in black or blue



Prod. No.		Size	Dimensions	
Black	Blue		Above patella	Below patella
00-3474-001-00	00-3474-006-00	Small	17 in. - 19 in. (43cm - 48cm)	12 1/2 in. - 14 in. (32cm - 36cm)
00-3474-002-00	00-3474-007-00	Medium	19 in. - 21 in. (48cm - 53cm)	14 in. - 15 1/2 in. (36cm - 39cm)
00-3474-003-00	00-3474-008-00	Large	21 in. - 23 in. (53cm - 58cm)	15 in. - 16 1/2 in. (39cm - 42cm)
00-3474-004-00	00-3474-009-00	X-Large	23 in. - 25 in. (58cm - 63cm)	16 in. - 17 1/2 in. (42cm - 44cm)
00-3474-005-00	00-3474-010-00	XX-Large	25 in. - 27 in. (63cm - 69cm)	17 in. - 18 1/2 in. (44cm - 47cm)

*Trademark of Royce Medical Company

Measurement is taken 6 in. (15cm) above/below patella

ZIMMER® PREMIUM NEOPRENE HINGED KNEE BRACE

- Dual action hinges incorporate the 3-dimensional molding technology
- Hinges are molded into the brace to maintain a low profile
- Accordion molding in the popliteal region helps eliminate bunching and provides exceptional patient comfort
- Lycra liner offers a breathable and soft contact area for comfortable and extensive use
- Varied areas of thickness give additional protection
- Available in black only



Prod. No.	Size	Knee Circumference
00-3475-010-00	Small	11 1/2 in. - 13 in. (29cm - 33cm)
00-3475-011-00	Medium	13 in. - 14 1/2 in. (33cm - 37cm)
00-3475-012-00	Large	14 1/2 in. - 16 in. (37cm - 41cm)
00-3475-013-00	X-Large	16 in. - 17 1/2 in. (41cm - 44cm)
00-3475-014-00	XX-Large	17 1/2 in. - 19 in. (44cm - 48cm)

Measured around the center of the knee cap with the leg fully extended

KNEE BRACES AND SUPPORTS

KNEE SLEEVE WITH CLOSED PATELLA AND PAD

- 1/8 in. (3mm) thick neoprene
- 13 in. (33cm) sleeve, anatomically contoured for better fit
- Provides warmth and moderate compression over the entire length of the support for minor knee injuries
- Available in beige or blue



Prod. No.			
Beige	Blue	Size	Knee Circumference
00-3455-003-00	00-3455-023-00	Small	13 in. - 14 in. (33cm - 36cm)
00-3455-004-00	00-3455-024-00	Medium	14 in. - 15 in. (36cm - 38cm)
00-3455-005-00	00-3455-025-00	Large	15 in. - 17 in. (38cm - 41cm)
00-3455-006-00	00-3455-026-00	X-Large	16 in. - 18 in. (41cm - 46cm)
00-3455-007-00	00-3455-027-00	XX-Large	18 in. - 20 in. (46cm - 51cm)

Measured around the center of the knee cap with the leg fully extended

KNEE SLEEVE WITH OPEN PATELLA

- 1/8 in. (3mm) thick neoprene
- Patellar cut-out relieves pressure over patella
- Quick and simple application
- Provides warmth and equal compression for minor knee injuries
- Available in beige or blue



Prod. No.			
Beige	Blue	Size	Knee Circumference
00-3456-003-00	00-3456-023-00	Small	13 in. - 14 in. (33cm - 36cm)
00-3456-004-00	00-3456-024-00	Medium	14 in. - 15 in. (36cm - 38cm)
00-3456-005-00	00-3456-025-00	Large	15 in. - 17 in. (38cm - 41cm)
00-3456-006-00	00-3456-026-00	X-Large	16 in. - 18 in. (41cm - 46cm)
00-3456-007-00	00-3456-027-00	XX-Large	18 in. - 20 in. (46cm - 51cm)

Measured around the center of the knee cap with the leg fully extended

UNIVERSAL KNEE SUPPORT WITH ADJUSTABLE BUTTRISS

- 1/8 in. (3mm) thick neoprene
- Support may be inverted for medial or lateral pull
- Removable buttress
- Patellar cut-out relieves pressure over patella
- Beige color



Prod. No.	Size
00-3457-000-00	Universal



KNEE BRACES AND SUPPORTS

PATELLAR KNEE SUPPORT WITH INFERIOR, SUPERIOR, RIGHT OR LEFT BUTTRESS

- Choice of 4 buttress positions to aid in treatment of patellar dislocations, tendonitis, subluxation and chondromalacia
- 1/8 in. (3mm) thick neoprene
- Patellar cut-out relieves pressure over patella
- Sewn-in horseshoe buttress made of 1/2 in. (13mm) diameter tubular neoprene borders patellar opening for stabilization



WITH INFERIOR BUTTRESS

Prod. No.		Size	Knee Circumference
Beige	Blue		
00-3458-103-00	00-3458-123-00	Small	13 in. - 14 in. (33cm - 36cm)
00-3458-104-00	00-3458-124-00	Medium	14 in. - 15 in. (36cm - 38cm)
00-3458-105-00	00-3458-125-00	Large	15 in. - 17 in. (38cm - 41cm)
00-3458-106-00	00-3458-126-00	X-Large	16 in. - 18 in. (41cm - 46cm)
00-3458-107-00	00-3458-127-00	XX-Large	18 in. - 20 in. (46cm - 51cm)

WITH SUPERIOR BUTTRESS

Prod. No.		Size	Knee Circumference
Beige	Blue		
00-3458-203-00	00-3458-223-00	Small	13 in. - 14 in. (33cm - 36cm)
00-3458-204-00	00-3458-224-00	Medium	14 in. - 15 in. (36cm - 38cm)
00-3458-205-00	00-3458-225-00	Large	15 in. - 17 in. (38cm - 41cm)
00-3458-206-00	00-3458-226-00	X-Large	16 in. - 18 in. (41cm - 46cm)
00-3458-207-00	00-3458-227-00	XX-Large	18 in. - 20 in. (46cm - 51cm)

WITH RIGHT BUTTRESS

Prod. No.		Size	Knee Circumference
Beige	Blue		
00-3458-303-00	00-3458-323-00	Small	13 in. - 14 in. (33cm - 36cm)
00-3458-304-00	00-3458-324-00	Medium	14 in. - 15 in. (36cm - 38cm)
00-3458-305-00	00-3458-325-00	Large	15 in. - 17 in. (38cm - 41cm)
00-3458-306-00	00-3458-326-00	X-Large	16 in. - 18 in. (41cm - 46cm)
00-3458-307-00	00-3458-327-00	XX-Large	18 in. - 20 in. (46cm - 51cm)

WITH LEFT BUTTRESS

Prod. No.		Size	Knee Circumference
Beige	Blue		
00-3458-403-00	00-3458-423-00	Small	13 in. - 14 in. (33cm - 36cm)
00-3458-404-00	00-3458-424-00	Medium	14 in. - 15 in. (36cm - 38cm)
00-3458-405-00	00-3458-425-00	Large	15 in. - 17 in. (38cm - 41cm)
00-3458-406-00	00-3458-426-00	X-Large	16 in. - 18 in. (41cm - 46cm)
00-3458-407-00	00-3458-427-00	XX-Large	18 in. - 20 in. (46cm - 51cm)

Measured around the center of the knee cap with the leg fully extended

POST-ARTHROSCOPY KNEE SUPPORT

- 1/8 in. (3mm) thick neoprene
- Adjustable patellar pads for added compression
- Provides compression and support
- Retains body warmth
- Beige color only



Prod. No. Size Knee Circumference

00-3461-003-00	Small	13 in. - 14 in. (33cm - 36cm)
00-3461-004-00	Medium	14 in. - 15 in. (36cm - 38cm)
00-3461-005-00	Large	15 in. - 17 in. (38cm - 41cm)
00-3461-006-00	X-Large	16 in. - 18 in. (41cm - 46cm)
00-3461-007-00	XX-Large	18 in. - 20 in. (46cm - 51cm)

Measured around the center of the knee cap with the leg fully extended

KNEE BRACES AND SUPPORTS

PATELLAR STABILIZED KNEE SUPPORT

- 1/8 in. (3mm) thick neoprene
- Patellar cut-out relieves pressure over patella
- Straps provide direct support to patellar area for added stability



Prod. No.			
Beige	Blue	Size	Knee Circumference
00-3459-003-00	00-3459-023-00	Small	13 in. - 14 in. (33cm - 36cm)
00-3459-004-00	00-3459-024-00	Medium	14 in. - 15 in. (36cm - 38cm)
00-3459-005-00	00-3459-025-00	Large	15 in. - 17 in. (38cm - 41cm)
00-3459-006-00	00-3459-026-00	X-Large	16 in. - 18 in. (41cm - 46cm)
00-3459-007-00	00-3459-027-00	XX-Large	18 in. - 20 in. (46cm - 51cm)

Measured around the center of the knee cap with the leg fully extended

HINGED KNEE SUPPORT

- Indicated for posterior cruciate instability. Also medial and lateral collateral ligament instability
- 1/8 in. (3mm) thick neoprene
- Exposed horseshoe superior buttress pad



Prod. No.			
Beige	Blue	Size	Knee Circumference
00-3460-003-00	00-3460-023-00	Small	13 in. - 14 in. (33cm - 36cm)
00-3460-004-00	00-3460-024-00	Medium	14 in. - 15 in. (36cm - 38cm)
00-3460-005-00	00-3460-025-00	Large	15 in. - 17 in. (38cm - 41cm)
00-3460-006-00	00-3460-026-00	X-Large	16 in. - 18 in. (41cm - 46cm)
00-3460-007-00	00-3460-027-00	XX-Large	18 in. - 20 in. (46cm - 51cm)

Measured around the center of the knee cap with the leg fully extended



LOWER LEG/ANKLE SUPPORTS

ANKLE SUPPORT WITH CONTACT CLOSURE

- General support for stiff or weakened ankles
- 1/8 in. (3mm) thick neoprene
- Heel cut-out eliminates potential skin irritation
- Easy to apply and remove
- Beige color



Prod. No.	Size	Shoe Size
00-3453-003-00	Small	6 - 7
00-3453-004-00	Medium	7 - 9
00-3453-005-00	Large	10 - 12
00-3453-006-00	X-Large	12 - 14

ZIMMER® GEL ANKLE BRACE

- D-ring lacing system helps add more compression and stability to help prevent inversion and eversion
- Anatomically designed shells lined with gel and foam pads help provide ankle protection and patient comfort
- Choice of gel pads available in full (80% gel/20% foam) or half lengths (50% gel/50% foam)
- Gel pads can be detached and placed in a freezer for cold therapy or warmed in the microwave for hot therapy
- Rigid plastic stays provide medial/lateral stability
- Pony size version is also available for children or small adults



GEL ANKLE BRACE 50/50		
Left	Right	Size
00-1421-001-00	00-1421-002-00	Adult, Universal
00-1423-001-00	00-1423-002-00	Pony



GEL ANKLE BRACE 80/20		
Left	Right	Size
00-1422-001-00	00-1422-002-00	Adult, Universal



ZIMMER® LACE-LOCK®* FOAM ANKLE BRACE

- Open cell pads help provide a lightweight, breathable brace
- D-ring lacing system helps add more compression and stability to help prevent inversion and eversion
- 3-dimensional molding technology varies padding for superior comfort and fit



Prod. No.		
Left	Right	Size
00-1424-001-00	00-1424-002-00	Adult, Universal
00-1425-001-00	00-1425-002-00	Pony

*Trademark of Royce Medical Company



LOWER LEG/ANKLE SUPPORTS

ZIMMER® PREMIUM HONEYCOMB ANKLE STIRRUP

- Built-in honeycomb comfort zones provide varying levels of compression and support
- Aeration holes ensure a cool and breathable pad surface
- Pivoting heel straps provide unobstructed range of motion in plantar and dorsiflexion
- Heel straps use contact closures to allow for quick adjustment and a custom fit even after the brace is applied
- Flex edge overmold technology incorporates a soft flexible edge onto the hard outer shell helping to eliminate pressure points on the ankle and increase comfort



Prod. No.		
Left	Right	Size
00-1435-001-00	00-1435-002-00	Universal

ZIMMER® AIR ANKLE STIRRUP

- Active air chambers deliver graduated pneumatic compression to reduce swelling and pain
- Pulsating multi-chamber bladder allows controlled displacement of air to provide a custom fit with each step
- Uniform compression provides increased stability and patient comfort
- Pivoting heel straps provide unobstructed range of motion in plantar and dorsiflexion
- Heel strap uses contact closures to allow for quick adjustment and a custom fit even after the brace is applied
- Flex edge overmold technology incorporates a soft flexible edge onto the hard outer shell helping to eliminate pressure points on the ankle and increase comfort
- Air bladders are pre-inflated and do not require pump inflating or straws reducing the risk of over-inflation



Prod. No.		
Left	Right	Size
00-1435-021-00	00-1435-022-00	Universal



LOWER LEG/ANKLE SUPPORTS

ZIMMER® EXO-ANKLE BRACE AND EXO-ANKLE BRACE WITH STRAPS

- Styles with or without nylon “figure 8” straps available
- Speed lacing system for easy application and adjustability of compression
- Breathable, lightweight fabric enhances comfort by allowing heat and moisture to escape the brace
- Slim design fits easily in athletic shoes
- Bi-lateral design fits left or right ankle
- Plastic exoskeleton provides stability and reduces ankle inversion
- Machine washable with air dry



WITHOUT FIGURE 8 STRAPS

Prod. No.	Size	Shoe Size	
		Women's	Men's
00-1436-001-00	Small	7 1/2 - 9	6 1/2 - 8
00-1436-002-00	Medium	9 1/2 - 11	8 1/2 - 12
00-1436-003-00	Large	11 1/2 - 13	12 1/2 - 14
00-1436-004-00	X-Large	13 1/2 - 15	

WITH FIGURE 8 STRAPS

Prod. No.	Size	Shoe Size	
		Women's	Men's
00-1436-005-00	Small	7 1/2 - 9	6 1/2 - 8
00-1436-006-00	Medium	9 1/2 - 11	8 1/2 - 12
00-1436-007-00	Large	11 1/2 - 13	12 1/2 - 14
00-1436-008-00	X-Large	13 1/2 - 15	

DELUXE ANKLE BRACE

- Contoured steel stays to fit malleolus and instep
- Perforated vinyl allows increased air circulation for patient comfort
- Contact closures for easy application and adjustment



Prod. No.	Size	Dimensions
00-1767-003-00	Small	7 in. - 8 1/2 in. (18cm - 22cm)
00-1767-004-00	Medium	8 1/2 in. - 10 in. (22cm - 25cm)
00-1767-005-00	Large	10 in. - 11 1/2 in. (25cm - 29cm)
00-1767-006-00	X-Large	11 1/2 in. - 13 in. (29cm - 33cm)

ANKLE BRACE WITH ELASTIC STRAP

- Heel-lock elastic strap provides increased medial/lateral support; restricts plantar/dorsal foot movement
- Rehabilitative support after grade I or grade II sprains



Prod. No.		Size	Shoe Size	
Left	Right		Women's	Men's
00-0766-001-00	00-0766-002-00	X-Small	Under 9	Under 8
00-0766-003-00	00-0766-004-00	Small	9 - 10	8 - 9
00-0766-005-00	00-0766-006-00	Medium	10 - 12	9 - 11
00-0766-007-00	00-0766-008-00	Large	12 - 14	11 - 13
00-0766-009-00	00-0766-010-00	X-Large	Over 14	Over 13

WALKERS

ZIMMER® FIXED HINGE LOWER LEG WALKER

- Contoured strut design allows the walker frame to conform to patient's anatomy ensuring a better fit
- Shock absorbing sole reduces impact from the heel strike to improve comfort during ambulation
- Rocker bottom for plantar unloading is engineered to be low profile and wide to minimize leg length discrepancy and help promote a natural gait with increased stability
- Plastic uprights provide support and increased protection
- Walker provides support and comfort during ambulation and replaces a walking cast
- Indicated for soft tissue injuries (grade II and III sprains), stable fractures, post-operative use, trauma, rehabilitation and other conditions recommended by a physician
- Lightweight design is easier on patients



WALKER

Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-004-00	Small	6 - 8	4 1/2 - 7
00-8916-005-00	Medium	8 1/2 - 11 1/2	7 1/2 - 10 1/2
00-8916-006-00	Large	11 1/2 - 13 1/2	10 1/2 - 12 1/2

REPLACEMENT LINERS

Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-004-01	Small	6 - 8	4 1/2 - 7
00-8916-005-01	Medium	8 1/2 - 11 1/2	7 1/2 - 10 1/2
00-8916-006-01	Large	11 1/2 - 13 1/2	10 1/2 - 12 1/2

ZIMMER® LOW TOP ANKLE WALKER

- Lowered and contoured strut for when a taller version is not needed
- Rocker bottom for plantar unloading is engineered to be low profile and wide to minimize leg length discrepancy and help promote a natural gait with increased stability



Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-041-00	Small	6 - 8	4 1/2 - 7
00-8916-042-00	Medium	8 1/2 - 11 1/2	7 1/2 - 10 1/2
00-8916-043-00	Large	11 1/2 - 13 1/2	10 1/2 - 12 1/2



WALKERS

ZIMMER® AIR WALKER

- Pneumatic system increases stability while decreasing pain and swelling
- Easy to use pump for inflation/deflation
- Air is easily adjusted to achieve desired compression and custom fit for better patient compliance
- Accommodates the swelling patterns that occur throughout the rehabilitation process as well as dressings
- Rocker bottom is specially engineered to be low and wide to help promote a natural gait and increased stability



Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-024-00	Small	6 - 8	4 1/2 - 7
00-8916-025-00	Medium	8 1/2 - 11 1/2	7 1/2 - 10 1/2
00-8916-026-00	Large	11 1/2 - 13 1/2	10 1/2 - 12 1/2

ZIMMER® LOWER LEG ORTHOSIS II

- Hinged ankle design provides up to 45 degrees of controlled ankle motion
- Low profile design provides better anatomical conformity for smoother gait
- Replacement liners are available



WALKER			
Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-014-00	Small	6 1/2 - 9 1/2	Up to 7
00-8916-015-00	Medium	10 - 12 1/2	8 - 9
00-8916-016-00	Large	13 - 15	10 and up

REPLACEMENT LINERS			
Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-011-00	Small	6 1/2 - 9 1/2	Up to 7
00-8916-012-00	Medium	10 - 12 1/2	8 - 9
00-8916-013-00	Large	13 - 15	10 and up

DH PRESSURE RELIEF™* WALKER

- Pre-inflated bladder for strut protection
- 90 degrees fixation of ankle to reduce shearing
- Selective pressure relief easily created by removal of independent shock-absorbing hexagons
- Rocker bottom is specially engineered to be low and wide to help promote a natural gait and increased stability
- Effective solution for treating mild to moderate plantar ulcers
- Inner bootie designed to consistently locate the foot relative to the pressure relief insole
- Patient-friendly leg straps
- Liner and sole are easily removed for cleaning



Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-031-00	Small	Up to 10	Up to 8
00-8916-032-00	Medium	10 1/2 - 12 1/2	8 1/2 - 10 1/2
00-8916-033-00	Large	12 1/2 and up	10 1/2 and up

*Trademark of Royce Medical Company

CAST SHOES/POST-OP SHOES

ROCKER SOLE CAST SHOE

- Curved rocker sole for easier walking
- Sole design promotes secure footing
- Durable denim upper



Prod. No.	Size	Dimensions
00-8716-011-00	X-Small	7 1/4 in. x 3 3/4 in. (18cm x 10cm)
00-8716-012-00	Small	9 1/4 in. x 4 1/4 in. (23cm x 11cm)
00-8716-013-00	Medium	10 7/8 in. x 5 in. (28cm x 13cm)
00-8716-014-00	Large	12 3/8 in. x 5 1/2 in. (31cm x 14cm)
00-8716-015-00	X-Large	13 1/2 in. x 5 7/8 in. (34cm x 15cm)

CUSHION ROCKER® CAST SHOE

- Rocker sole promotes normal gait
- Tread designed for secure footing
- Contact closure straps for easy application and removal



Prod. No.	Size	Dimensions
00-1995-009-00	Pediatric	2 1/2 in. x 4 1/2 in. (6cm x 11cm)
00-1995-010-00	XX-Small	3 in. x 5 in. (8cm x 13cm)
00-1995-011-00	X-Small	7 1/4 in. x 3 3/4 in. (18cm x 10cm)
00-1995-012-00	Small	9 1/4 in. x 4 1/4 in. (23cm x 11cm)
00-1995-013-00	Medium	10 7/8 in. x 5 in. (28cm x 13cm)
00-1995-014-00	Large	12 3/8 in. x 5 1/2 in. (31cm x 14cm)
00-1995-015-00	X-Large	13 1/2 in. x 5 7/8 in. (34cm x 15cm)

ZIMMER® CUSHION ROCKER® CAST BOOT

- Closed toe design for complete protection against dirt and moisture
- Rocker sole promotes normal gait
- Durable canvas material construction



Prod. No.	Size	Dimensions
00-1995-021-00	X-Small	7 1/4 in. x 3 3/4 in. (18cm x 10cm)
00-1995-022-00	Small	9 1/4 in. x 4 1/4 in. (23cm x 11cm)
00-1995-023-00	Medium	10 7/8 in. x 5 in. (28cm x 13cm)
00-1995-024-00	Large	12 3/8 in. x 5 1/2 in. (31cm x 14cm)
00-1995-025-00	X-Large	13 1/2 in. x 5 7/8 in. (34cm x 15cm)

DARCO®* CAST BOOT

- Neoprene material for longer wear and snug fit
- Contact closures simplify application and adjustment
- Rocker sole promotes patient mobility
- Available in navy blue



Prod. No.	Size	Dimensions
00-1997-001-00	X-Small	8.8 in. x 3.1 in. (22cm x 8cm)
00-1997-002-00	Small	9.4 in. x 3.3 in. (24cm x 8cm)
00-1997-003-00	Medium	10.4 in. x 3.7 in. (26cm x 9cm)
00-1997-004-00	Large	11.2 in. x 3.9 in. (28cm x 10cm)
00-1997-005-00	X-Large	12 in. x 4.1 in. (30cm x 10cm)

*Trademark of Darco International, Inc.



CAST SHOES/POST-OP SHOES

ZIMMER® POST-OP SHOE

- Accommodates most digital, forefoot, full foot, or ankle dressings
- Perforated vinyl allows increased air circulation for patient comfort
- Completely padded for maximum comfort
- Slip resistant outer sole is comprised of corrugated rubber bonded to wood
- Insole is foam
- Lace closures ensure secure fit



WOMEN'S

Prod. No.	Size	Dimensions
00-1991-004-00	Small	3 1/2 in. x 9 1/2 in. (9cm x 24cm)
00-1991-005-00	Medium	3 3/4 in. x 10 in. (10cm x 25cm)
00-1991-006-00	Large	3 7/8 in. x 10 1/2 in. (10cm x 27cm)

MEN'S

Prod. No.	Size	Dimensions
00-1991-001-00	Small	3 7/8 in. x 11 1/4 in. (10cm x 29cm)
00-1991-002-00	Medium	4 1/8 in. x 11 1/2 in. (10cm x 29cm)
00-1991-003-00	Large	4 3/8 in. x 12 in. (11cm x 30cm)

PEDIATRIC

Prod. No.	Size	Dimensions
00-8601-697-01	Pediatric	2 1/2 in. x 5 in. (6cm x 13cm)
00-8601-697-02	Youth	3 in. x 8 in. (8cm x 20cm)

ZIMMER® PADDED POST-OP SHOE

- Inner sole accommodates natural arch of foot
- Padded tongue reduces pressure on dressing and sore areas
- Perforated Tietex material allows increased air circulation for patient comfort
- Convenient contact closures help allow easy application



WOMEN'S

Prod. No.	Size	Dimensions
00-1611-024-00	Small	3 1/16 in. x 9 3/8 in. (8cm x 24cm)
00-1611-025-00	Medium	3 1/4 in. x 10 in. (8cm x 25cm)
00-1611-026-00	Large	4 3/16 in. x 12 3/16 in. (11cm x 31cm)

MEN'S

Prod. No.	Size	Dimensions
00-1611-021-00	Small	4 in. x 10 3/4 in. (10cm x 27cm)
00-1611-022-00	Medium	4 1/16 in. x 11 3/8 in. (10cm x 29cm)
00-1611-023-00	Large	4 3/16 in. x 12 3/16 in. (11cm x 31cm)

DARCO®* POST-OP SHOE

- Forefoot support that provides more rigid control and better protection after osteotomies
- Sole angled at 15 degrees declination for reduced heel strike
- Reinforced heel cup helps stabilize foot by reducing pronation
- Polywood mid-sole provides semi-rigidity and control
- Breathable mesh upper conforms to bony abnormalities
- Padded heel collar helps to reduce slippage
- Extra-long straps accommodate even the bulkiest dressings



WOMEN'S

Prod. No.	Size	Women's Shoe Size
00-1711-004-00	Small	4 - 6
00-1711-005-00	Medium	6 1/2 - 8
00-1711-006-00	Large	8 1/2 - 10

MEN'S

Prod. No.	Size	Men's Shoe Size
00-1711-011-00	Small	6 - 8
00-1711-012-00	Medium	8 1/2 - 10
00-1711-013-00	Large	10 1/2 - 12
00-1711-014-00	X-Large	12 1/2 - 14

*Trademark of Darco International, Inc.

CAST SHOES/POST-OP SHOES

DARCO®* MED/SURG POST-OP SHOE

- Padded fold-over contact closure for patient comfort and easy application and removal
- Semi-rigid sole with greater support and better protection after surgery or injury
- Breathable mesh upper conforms to bony abnormalities of the foot
- Reinforced heel counter and padded heel collar to stabilize the rear foot and reduce slippage
- Increased tongue area allows for large bulky dressings
- Classic black design



Prod. No.	Size	Shoe Size	
		Women's	Men's
00-1811-001-00	X-Small	3 - 6 1/2	
00-1811-002-00	Small	7 - 9 1/2	5 - 6 1/2
00-1811-003-00	Medium		7 - 8 1/2
00-1811-004-00	Large		9 - 11
00-1811-005-00	X-Large		11 1/2 - 14

*Trademark of Darco International, Inc.

DARCO®* ORTHOWEDGE™* HEALING SHOE

- Strapless forefoot closure eliminates pressure point caused by buckles
- Short wedge offloads an even larger area of the forefoot for greater pressure relief
- Heel locking strap eliminates slippage even after partial amputations
- Weight is redistributed behind metatarsal heads and digits during the weight-bearing portion of the gait cycle
- Indicated for the treatment of diabetic ulcerations and post-op after surgeries to correct hallux valgus, hammer toes, Taylor's bunions and other osteotomies with or without pins
- Classic black design



Prod. No.	Size	Shoe Size	
		Women's	Men's
00-1443-001-00	X-Small	4 - 7	
00-1443-002-00	Small	7 1/2 - 10	6 - 8
00-1443-003-00	Medium		8 1/2 - 10
00-1443-004-00	Large		10 1/2 - 12
00-1443-005-00	X-Large		12 1/2 - 14

*Trademarks of Darco International, Inc.



CAST SHOES/POST-OP SHOES

DH PRESSURE RELIEF™* SHOE

- Selective pressure relief easily created by removal of independent shock-absorbing hexagons
- Discreet black shoe helps to ensure compliance
- Removable toe piece allows for a sandal-look or shoe-look
- Sole is easily removed for cleaning
- Rocker bottom helps to promote a natural gait and increased stability



Prod. No.	Size	Shoe Size	
		Women's	Men's
00-1433-001-00	Small	6 1/2 - 8 1/2	4 - 6
00-1433-002-00	Medium	9 1/2 - 11	6 1/2 - 8 1/2
00-1433-003-00	Large	12 and up	9 1/2 - 11
00-1433-004-00	X-Large		11 1/2 - 14

*Trademark of Royce Medical Company

MISCELLANEOUS SUPPORTS

PAVLIK HARNESS

- Provides support for hip displasia and dislocation without discomfort of cast
- Allows necessary joint movement while child achieves self-reduction
- Washable
- Contact closures for easy application and removal



Prod. No.	Size	Chest Band Range
00-8601-416-00	X-Small	10 1/2 in. x 14 in. (27cm x 36cm)
00-1716-010-00	Small	13 1/2 in. x 16 1/2 in. (34cm x 42cm)
00-1716-011-00	Medium	15 1/4 in. x 19 1/4 in. (39cm x 49cm)
00-1716-012-00	Large	18 3/4 in. x 22 3/4 in. (48cm x 58cm)

ZIMMER® PREMIUM NIGHT SPLINT

- Ideal for plantar fasciitis, achilles tendonitis and other lower extremity overuse injuries
- 3-dimensional molded technology allows for variable padding thickness in strategic locations for superior comfort and fit
- Soft, lycra-lined foam interior is plush and breathable helping to eliminate heat and moisture build-up
- Lightweight, low profile shell is sturdy
- Padded straps with easy-release buckles ensure immobilization while maintaining comfort across the flexor tendons
- Adjustable dorsiflexion straps provide a gentle stretching of the plantar fascia and achilles tendon
- Adjustable from 90 degrees to 10 degrees dorsiflexion



Prod. No.	Size	Shoe Size	
		Women's	Men's
00-1444-001-00	Small	Up to 7 1/2	Up to 7
00-1444-002-00	Medium	8 - 10 1/2	7 1/2 - 10
00-1444-003-00	Large	11 and up	10 1/2 and up

ORTHOPAEDIC WRAPS

ZIMMER® ORTHOPAEDIC WRAPS

- Cotton/polyester knit construction for conformity, durability and comfort
- Lightweight, breathable material
- Latex dipped ends (not sewn) won't unravel
- Premium wrap available sterile or non-sterile
- Non-latex wrap also available
- 4 separate widths
- Firm, consistent support and compression to reduce swelling and relieve pain associated with sprains and strains
- Washable (retains elasticity after repeated use and washings)
- 2 in. (5cm) for wrist and feet, 3 in. (8cm) for elbow, ankle, knee and lower leg, 4 in. (10cm) for knee, lower leg and shoulder, 6 in. (15cm) shoulder and chest



ZIMMER PREMIUM ORTHOPAEDIC WRAP

Latex	Non-Latex	Width/Length	Packaging
00-5555-002-00	00-5555-102-00	2 in. x 5.5 yds. (5cm x 5m)	10/box
00-5555-003-00	00-5555-103-00	3 in. x 5.5 yds. (8cm x 5m)	10/box
00-5555-004-00	00-5555-104-00	4 in. x 5.5 yds. (10cm x 5m)	10/box
00-5555-006-00	00-5555-106-00	6 in. x 5.5 yds. (15cm x 5m)	10/box
00-5555-024-00	00-5555-124-00	4 in. x 11 yds. (10cm x 10m)	10/box
00-5555-026-00	00-5555-126-00	6 in. x 11 yds. (15cm x 10m)	10/box

ZIMMER STERILE PREMIUM ORTHOPAEDIC WRAP

Latex	Non-Latex	Width/Length	Packaging
00-5555-032-00	00-5555-232-00	2 in. x 5.5 yds. (5cm x 5m)	10/box
00-5555-033-00	00-5555-233-00	3 in. x 5.5 yds. (8cm x 5m)	10/box
00-5555-034-00	00-5555-234-00	4 in. x 5.5 yds. (10cm x 5m)	10/box
00-5555-036-00	00-5555-236-00	6 in. x 5.5 yds. (15cm x 5m)	10/box
00-5555-044-00	00-5555-244-00	4 in. x 11 yds. (10cm x 10m)	10/box
00-5555-046-00	00-5555-246-00	6 in. x 11 yds. (15cm x 10m)	10/box

ZIMMER STANDARD ORTHOPAEDIC WRAP

Latex	Non-Latex	Width/Length	Packaging
00-4444-002-00	00-4444-102-00	2 in. x 5 yds. (5cm x 5m)	10/box
00-4444-003-00	00-4444-103-00	3 in. x 5 yds. (8cm x 5m)	10/box
00-4444-004-00	00-4444-104-00	4 in. x 5 yds. (10cm x 5m)	10/box
00-4444-006-00	00-4444-106-00	6 in. x 5 yds. (15cm x 5m)	10/box
00-4444-012-00	00-4444-112-00	2 in. x 5 yds. (5cm x 5m)	100/box
00-4444-013-00	00-4444-113-00	3 in. x 5 yds. (8cm x 5m)	100/box
00-4444-014-00	00-4444-114-00	4 in. x 5 yds. (10cm x 5m)	100/box
00-4444-016-00	00-4444-116-00	6 in. x 5 yds. (15cm x 5m)	50/box
00-4444-024-00	00-4444-124-00	4 in. x 10 yds. (10cm x 9m)	10/box
00-4444-026-00	00-4444-126-00	6 in. x 10 yds. (15cm x 9m)	10/box



COLD THERAPY

ZIMMER® COLD THERAPY WRAPS

- Place reusable cold packs in freezer for 2 to 4 hours to activate
- Insert the pack in the neoprene wrap and apply to injury site with an article of clothing or towel between the skin and the neoprene wrap for comfortable cold therapy
- Cold packs remain flexible to better conform to the patient's anatomy and injury site for maximum therapeutic benefit
- Stretchable wrap is secured with a contact closure allowing for varying degrees of compression and easy application
- Combination of cold and compression quickly constrict blood vessels to control pain and swelling



Prod. No.	Size	Description	Dimensions
00-2761-022-12	X-Small	Wrist and Forearm	22 in. x 3 1/2 in. (56cm x 9cm)
00-2761-023-12	Small	Ankle and Elbow	23 in. x 6 1/2 in. (58cm x 17cm)
00-2761-024-12	Medium	Ankle, Elbow and Knee	31 in. x 7 in. (79cm x 18cm)
00-2761-025-12	Large	Knee, Thigh and Back	40 in. x 8 1/2 in. (102cm x 22cm)
00-2761-026-12	X-Large	Thigh and Back	47 in. x 11 1/2 in. (119cm x 29cm)
00-2763-025-12	X-Long	Thigh and Back	47 in. x 8 1/2 in. (119cm x 22cm)
00-2761-027-12	Universal	Shoulder	23 1/4 in. x 11 in. (59cm x 28cm)

ZIMMER® INSTANT COLD PACK

- Indicated for cold therapy to reduce pain and/or swelling
- Squeezing pouch produces instant cold (32°F) that lasts for approximately 30 minutes
- Remains soft and pliable to conform to body contours



Prod. No.	Dimensions	Packaging
00-1610-014-00	4 3/8 in. x 7 1/8 in. (11cm x 18cm)	32 pks/box
00-1610-015-00	5 3/4 in. x 9 1/8 in. (15cm x 23cm)	16 pks/box

PRESSURE CONTROLLED STOCKINGS

ZIMMER® PCS PRESSURE CONTROLLED STOCKINGS

- Clear view inspection window allows constant monitoring of toes for signs of circulation problems
- Looser knit at heel and instep helps avoid constriction and enhances patient comfort
- Color-coded packaging for easy size selection



KNEE-HIGH STOCKINGS—10 PAIR PER BOX

Prod. No.	Size	Calf Circumference	Length
00-2040-003-00	S, Short	Up to 12 in. (30cm)	Up to 14 in. (36cm)
00-2040-007-00	S, Regular		Up to 16 in. (41cm)
00-2040-011-00	S, Long		Over 16 in. (41cm)
00-2040-004-00	M, Short	12 in. - 15 in. (30cm - 38cm)	Up to 15 in. (38cm)
00-2040-008-00	M, Regular		Up to 17 in. (43cm)
00-2040-012-00	M, Long		Over 17 in. (43cm)
00-2040-005-00	L, Short	15 in. - 18 in. (38cm - 46cm)	Up to 16 in. (41cm)
00-2040-009-00	L, Regular		Up to 18 in. (46cm)
00-2040-013-00	L, Long		Over 18 in. (46cm)
00-2040-010-00	XL, Regular	Over 18 in. (46cm)	Up to 18 in. (46cm)
00-2040-014-00	XL, Long		Over 18 in. (46cm)

THIGH-HIGH STOCKINGS—5 PAIR PER BOX

Prod. No.	Size	Calf Circumference	Length
00-2040-023-00	S, Short	Up to 12 in. (30cm)	Up to 29 in. (74cm)
00-2040-027-00	S, Regular		Up to 33 in. (84cm)
00-2040-031-00	S, Long		Over 33 in. (84cm)
00-2040-024-00	M, Short	12 in. - 15 in. (30cm - 38cm)	Up to 29 in. (74cm)
00-2040-028-00	M, Regular		Up to 33 in. (84cm)
00-2040-032-00	M, Long		Over 33 in. (84cm)
00-2040-025-00	L, Short	15 in. - 18 in. (38cm - 46cm)	Up to 29 in. (74cm)
00-2040-029-00	L, Regular		Up to 33 in. (84cm)
00-2040-033-00	L, Long		Over 33 in. (84cm)
00-2040-030-00	XL, Regular	Over 18 in. (46cm)	Up to 33 in. (84cm)
00-2040-034-00	XL, Long		Over 33 in. (84cm)



SPLINTS

ZIMMER® FIBERGLASS SPLINTING SYSTEM

- Provides rigid external immobilization for fractures, sprains, dislocations and other skeletal deformities following surgical or non-surgical correction
- Features 6 layers of fiberglass tape covered with clear plastic film and soft foam padding
- Water activated



Prod. No.	Dimensions	Box Quantity
00-8906-013-00	3 in. x 12 in. (8cm x 30cm)	10
00-8906-014-00	4 in. x 15 in. (10cm x 38cm)	10
00-8906-024-00	4 in. x 30 in. (10cm x 76cm)	5

UNIVERSAL PLASTER SPLINTING SYSTEM

- 13 layers of extra fast setting plaster covered with open cell foam and soft cotton stockinette
- Time-saving—just trim to fit, immerse in water, mold to extremity and cover with elastic wrap. Sets in approximately 5 minutes



Prod. No.	Dimensions	Quantity
00-8905-002-00	2 in. x 260 in. (5cm x 66m)	1 roll
00-8905-003-00	3 in. x 260 in. (8cm x 66m)	1 roll
00-8905-004-00	4 in. x 260 in. (10cm x 66m)	1 roll
00-8905-005-00	5 in. x 260 in. (13cm x 66m)	1 roll
00-8905-006-00	6 in. x 260 in. (15cm x 66m)	1 roll

STAX FINGER SPLINTS

- Supports distal joint of finger in extension, yet permits unrestricted movement of proximal interphalangeal joint
- No padding required
- Available in 8 sizes



Prod. No.	Dimensions	Box Quantity
00-2549-001-00 (Kit)	Assortment of all sizes 1 in. - 7 in. (3cm - 18cm)	30
00-2549-011-00	1 in. (3cm)	4
00-2549-012-00	2 in. (5cm)	4
00-2549-013-00	3 in. (8cm)	4
00-2549-014-00	4 in. (10cm)	5
00-2549-015-00	5 in. (13cm)	4
00-2549-016-00	5 1/2 in. (14cm)	4
00-2549-017-00	6 in. (15cm)	4
00-2549-018-00	7 in. (18cm)	1

Replacement splints may be ordered in packages of 5.

ZIMFOAM® FINGER SPLINT KIT

- Assortment of 58 finger splints in convenient storage container
- Ideal for physician's offices, small clinics



Prod. No.	Description	Box Quantity
00-1756-001-00	Complete Kit. 58 assorted sizes	
00-0756-001-00	2 3/4 in. (7cm) Baseball Splint	6
00-1756-002-00	2 13/16 in. (7cm) Baseball Splint	6
66-0017-001-00	1 5/8 in. (4cm) Finger Splint—Padded	6
66-0017-002-00	2 5/8 in. (7cm) Finger Splint—Padded	6
66-0017-003-00	6 in. (15cm) Finger Splint—Padded	12
66-0017-004-00	1 1/4 in. (3cm) Finger Cot—Padded	6
66-0017-005-00	2 in. (5cm) Finger Cot—Padded	8
66-0017-006-00	3 1/2 in. (8cm) Finger Cot—Padded	8

Replacement splints may be ordered in boxes of 12.

SPLINTS

DIGITAL FINGER SPLINT KIT

- Aluminum splints easily shaped to flex or hyperextend joints
- Longer splints include a large circular base for stable splinting over metacarpal-phalangeal joint
- Assortment of 54 splints, 6 each of 9 sizes



Prod. No.	Description
00-2548-001-00	Complete kit. 54 assorted sizes.
00-2548-011-00	1 1/4 in. (3cm) Aluminum Finger Splint
00-2548-012-00	2 1/4 in. (6cm) Aluminum Finger Splint
00-2548-013-00	3 1/4 in. (8cm) Aluminum Finger Splint
00-2548-014-00	4 1/4 in. (11cm) Aluminum Finger Splint
00-2548-015-00	3 1/4 in. (8cm) Aluminum Finger Splint with Bulb
00-2548-016-00	4 1/4 in. (11cm) Aluminum Finger Splint with Bulb
00-2548-017-00	5 1/4 in. (13cm) Aluminum Finger Splint with Bulb
00-2548-018-00	6 1/4 in. (16cm) Aluminum Finger Splint with Bulb
00-2548-019-00	7 1/4 in. (18cm) Aluminum Finger Splint with Bulb

Replacement splints sold in packages of 6.

ZIMFOAM® BASEBALL SPLINT

- Maintains interphalangeal joints in correct position following distal joint injuries
- Malleable arms fold to hold finger—no adhesive required
- 1/4 in. (6mm) *Zimfoam* padding



Prod.No.	Length	Box Quantity
00-0756-011-00 Small	2 1/4 in. (6cm)	12
00-0756-001-00 Medium	2 3/4 in. (7cm)	12
00-0756-002-00 Large	2 13/16 in. (7cm)	12

ALUMINUM FINGER SPLINT

- Indicated for transverse fractures of the phalanges
- Curved to conform to finger in semi-flexed position
- Narrow width for women's or children's fingers
- Packaged in boxes of 12



Padded Regular Width	Padded Narrow Width	Length
66-0017-001-00	66-0017-011-00	1 5/8 in. (4cm)
66-0017-002-00	66-0017-012-00	2 5/8 in. (7cm)
66-0017-003-00	66-0017-013-00	6 in. (15cm)
66-0017-008-00		3 5/8 in. (9cm)
66-0017-009-00		4 5/8 in. (12cm)

FINGER COT

- Lightweight aluminum protects injuries to fingertips
- Provides good ventilation
- Both padded and non-padded styles



Padded	Length	Box Quantity
66-0017-004-00	1 1/4 in. (3cm)	2
66-0017-005-00	2 in. (5cm)	2
66-0017-006-00	3 in. (8cm)	2
66-0017-007-00	3 1/2 in. (9cm)	2



SPLINTS

ZIMFOAM® FINGER SPLINT


- Can be applied as finger cot, baseball splint or finger protector
- Easily formed or cut to desired shape with ordinary scissors
- Aluminum with 1/4 in. (6mm) *Zimfoam* padding



Prod. No.	Dimensions	Box Quantity
00-0548-001-00	18 in. x 3/4 in. (46cm x 2cm)	12
00-0548-002-00	18 in. x 1 in. (46cm x 3cm)	12
00-0548-005-00	9 in. x 1/2 in. (23cm x 1cm)	12

FINGER PROTECTOR

- Lightweight aluminum for patient comfort, convenience
- Available with or without padding



Non-Padded	Padded	Length	Box Quantity
00-0025-001-00	66-0025-001-00	1 1/2 in. (4cm)	12
00-0025-002-00	66-0025-002-00	3 in. (8cm)	12
00-0025-003-00	66-0025-003-00	3 1/2 in. (9cm)	12

LEWIN BASEBALL SPLINT

- Indicated for distal fractures of the distal phalanx
- Aluminum with 1/4 in. (6mm) padding



Prod. No.	Length	Box Quantity
66-0081-000-00	3 in. (8cm)	12
66-0081-002-00	4 in. (10cm)	12
66-0081-003-00	5 in. (13cm)	12

LEWIN STERN THUMB SPLINT

- Padded aluminum can be shaped to maintain thumb in flexion or extension
- Adjustable leather wrist band one size fits all



Prod. No.	Dimensions
66-0080-001-00	Wrist Band 6 3/8 in. (16cm) Splint 6 in. (15cm)

LEWIN STERN FINGER SPLINT

- Aluminum splint with *Zimfoam* padding
- Permits splinting of any finger in varying degrees of flexion
- Adjustable leather wrist band



Prod. No.	Dimensions
66-0080-002-00	Wrist Band 6 3/8 in. (16cm) Splint 8 7/16 in. (21cm)

SPLINTS

LYTLE METACARPAL SPLINT

- Ideal for metacarpal and phalangeal fractures
- Smooth anodized aluminum conforms to palm
- Also maintains slight amount of traction



Prod. No.	Size	Anterior	Length
00-0016-001-00	Youth, left	3 7/16 in. (9cm)	3 3/8 in. (9cm)
00-0016-002-00	Youth, right	3 7/16 in. (9cm)	3 3/8 in. (9cm)
00-0016-003-00	Adult, left	3 7/8 in. (10cm)	3 3/8 in. (9cm)
00-0016-004-00	Adult, right	3 7/8 in. (10cm)	3 3/8 in. (9cm)

MASON-ALLEN HAND SPLINT

- Immobilizes either hand in position of function following fracture or infection
- Helps prevent contractures in burn cases
- Conforms to convexity of forearm. Molded cup fits palmer concavity and heel of palm



Prod. No.	Size	Length
00-0376-001-00	Medium	10 7/16 in. (27cm)
00-0376-002-00	Large	12 1/4 in. (31cm)

ZIMMER® ALUMINUM COLLES SPLINT

- Perforated for air circulation
- Radiolucent
- *Zimfoam* padding included



Left	Right	Size	Dimensions
66-0020-001-00	66-0020-002-00	Small child	5 5/6 in. x 2 3/8 in. (15cm x 6cm)
66-0020-003-00	66-0020-004-00	Child	8 1/4 in. x 2 3/4 in. (21cm x 7cm)
66-0020-005-00	66-0020-006-00	Youth small adult	10 in. x 3 1/4 in. (25cm x 8cm)
66-0020-007-00	66-0020-008-00	Adult	11 in. x 4 in. (28cm x 10cm)

PLASTIC COLLES SPLINT

- Lightweight and economical
- Radiolucent
- *Zimfoam* padding included



Left	Right	Size	Dimensions
00-1704-005-00	00-1704-006-00	Small child	5 5/6 in. x 2 3/8 in. (15cm x 6cm)
00-1704-007-00	00-1704-008-00	Child	8 1/4 in. x 2 3/4 in. (21cm x 7cm)
00-1704-009-00	00-1704-010-00	Youth small adult	10 in. x 3 1/4 in. (25cm x 8cm)
00-1704-011-00	00-1704-012-00	Adult	11 in. x 4 in. (28cm x 10cm)

ALUMINUM FOREARM SPLINT

- For fractures of the radius and ulna
- Perforated for air circulation
- *Zimfoam* padding for comfort
- Radiolucent



Left	Right	Size	Length
66-0012-001-00	66-0012-002-00	Child	7 1/2 in. (19cm)
66-0012-003-00	66-0012-004-00	Youth	10 in. (25cm)
66-0012-005-00	66-0012-006-00	Adult	11 1/2 in. (29cm)



SPLINTS

ALUMINUM POSTERIOR TIBIA AND FIBULA SPLINTS

- Epoxy-coated aluminum with 1/4 in. (6mm) *Zimfoam* padding
- Available in hip and knee lengths with airholes for ventilation
- Radiolucent



KNEE LENGTH

Prod. No.	Size	Length
66-0040-001-00	Small child	8 3/4 in. (22cm)
66-0040-002-00	Child	11 3/4 in. (30cm)
66-0040-003-00	Youth small adult	14 in. (36cm)
66-0040-004-00	Adult	15 3/4 in. (40cm)

HIP LENGTH

Prod. No.	Size	Length
66-0040-005-00	Small child	14 in. (36cm)
66-0040-006-00	Child	19 1/2 in. (50cm)
66-0040-007-00	Youth small adult	25 1/4 in. (64cm)
66-0040-008-00	Adult	31 in. (79 cm)

PLASTIC POSTERIOR TIBIA AND FIBULA SPLINTS

- Available in choice of either knee or hip length sizes
- *Zimfoam* padding for comfort
- Radiolucent



KNEE LENGTH

Prod. No.	Size	Length
00-1706-005-00	Small child	8 3/4 in. (22cm)
00-1706-006-00	Child	11 5/8 in. (30cm)
00-1706-007-00	Youth	14 in. (36cm)
00-1706-008-00	Adult	15 3/4 in. (40cm)

HIP LENGTH

Prod. No.	Size	Length
00-1706-009-00	Small child	14 in. (36cm)
00-1706-010-00	Child	19 1/2 in. (50cm)
00-1706-011-00	Youth	25 1/4 in. (64cm)
00-1706-012-00	Adult	31 in. (79cm)

HUMERAL FUNCTIONAL FRACTURE BRACE

- Foam padded, perforated plastic sleeve with plastic tongue, nylon straps with contact closure. More comfortable and convenient than plaster casts
- Pre-formed braces eliminate time-consuming pattern cutting and molding
- Can be removed for extremity washing if approved by physician



Prod. No.	Size	Limb Circumference
00-1632-010-00	Small	7 in. - 11 in. (18cm - 28cm)
00-1632-011-00	Large	9 in. - 14 in. (23cm - 36cm)

SPLINTS

THUMB-WRIST IMMOBILIZER

- For stabilizing ulnar collateral ligament injuries, Colles fractures or sprained wrists
- Holds thumb in extended abducted position and wrist in neutral position
- Lightweight polyethylene with foam lining and contact closures



Prod. No.	Size	Width at MP Joint
00-8915-043-00 (Left)	Small	2 3/4 in. - 3 1/2 in. (7cm - 9cm)
00-8915-045-00 (Left)	Medium	3 1/2 in. - 4 in. (9cm - 10cm)
00-8915-047-00 (Left)	Large	4 in. - 4 1/2 in. (10cm - 11cm)
00-8915-044-00 (Right)	Small	2 3/4 in. - 3 1/2 in. (7cm - 9cm)
00-8915-046-00 (Right)	Medium	3 1/2 in. - 4 in. (9cm - 10cm)
00-8915-048-00 (Right)	Large	4 in. - 4 1/2 in. (10cm - 11cm)

CASTING PRODUCTS

ZIMMER® FIBERGLASS CASTING TAPE

- Available in four widths—2 in. - 5 in. (5cm - 13cm)
- Choice of 6 colors
- Cast gel included with each box for those who prefer a slippery tape during application
- Packaged 10 rolls per box
- Each roll is 4 yds. (4m) long
- Radiolucent



Prod. No.	Color	Width
00-8897-002-00	Regular	2 in. (5cm)
00-8897-002-01	Red	2 in. (5cm)
00-8897-002-02	Pink	2 in. (5cm)
00-8897-002-03	Light Blue	2 in. (5cm)
00-8897-002-04	Blue	2 in. (5cm)
00-8897-002-06	Green	2 in. (5cm)
00-8897-003-00	Regular	3 in. (8cm)
00-8897-003-01	Red	3 in. (8cm)
00-8897-003-02	Pink	3 in. (8cm)
00-8897-003-03	Light Blue	3 in. (8cm)
00-8897-003-04	Blue	3 in. (8cm)
00-8897-003-06	Green	3 in. (8cm)
00-8897-004-00	Regular	4 in. (10cm)
00-8897-004-01	Red	4 in. (10cm)
00-8897-004-02	Pink	4 in. (10cm)
00-8897-004-03	Light Blue	4 in. (10cm)
00-8897-004-04	Blue	4 in. (10cm)
00-8897-004-06	Green	4 in. (10cm)
00-8897-005-00	Regular	5 in. (13cm)



CASTING PRODUCTS

FIBERGLASS CAST SPLINTS

- Reinforces stress areas on a cast
- Helps ensure against breakdown
- Packaged 10 rolls per box



Product No.	Width	Length
00-8896-013-00	3 in. (8cm)	15 in. (38cm)
00-8896-014-00	4 in. (10cm)	15 in. (38cm)
00-8896-023-00	3 in. (8cm)	24 in. (61cm)
00-8896-024-00	4 in. (10cm)	24 in. (61cm)

SYNTHETIC STOCKINETTE

- 4 widths in handy dispenser cartons
- Polyester material
- Sold in 25 yds. (23m) lengths and choice of 4 widths
- Can be used with synthetic or plaster casts



Prod. No.	Width
00-8889-002-00	2 in. (5cm)
00-8889-003-00	3 in. (8cm)
00-8889-004-00	4 in. (10cm)
00-8889-005-00	5 in. (13cm)

SYNTHETIC CAST PADDING


- Cushions bony prominences and stress areas
- Self-adhering to simplify application
- Each roll is 4 yds. (4m) long, 20 rolls per box
- Choice of 4 widths



Prod. No.	Width
00-8890-002-00	2 in. (5cm)
00-8890-003-00	3 in. (8cm)
00-8890-004-00	4 in. (10cm)
00-8890-005-00	5 in. (13cm)

WHITE FELT PADDING

- Polyester padding for casts, splints, braces, prostheses
- Choice of 4 thicknesses
- Sold in rolls 36 in. wide by 21 in. long (91cm x 53cm)



Product No.	Thickness
00-0211-001-00	3/8 in. (10mm)
00-0217-005-00	1/8 in. (3mm)
00-0775-001-00	1/4 in. (6mm)
00-0775-002-00	1/2 in. (13mm)

ZIMMER® CAST CUTTER II

- High RPM, high-torque motor provides power for cutting all types of casting materials
- Handle designed for positive control and maneuverability
- Includes one 2 1/2 in. (6 cm) diameter Teflon coated universal blade and blade wrench



Prod. No.	Description
00-8909-001-00	Zimmer Cast Cutter II

CASTING PRODUCTS

CAST CUTTER BLADES

- Universal blade with Teflon coating for cutting both plaster and synthetic casts
- Stainless steel Premium blade features gold titanium nitride plating for longer wear
- Both blades 2 1/2 in. (6cm) in diameter



Product No.	Description	Packaging
00-8909-011-00	Premium Cast Cutter Blade	3/box
00-8909-012-00	Universal Cast Cutter Blade	6/box

ZIMMER® CAST REMOVAL SYSTEM

- Includes cast cutter, hose, vacuum and support frame
- Keeps cast room clean and dust free
- Filter system helps keep even the smallest particles out of motor and air



Prod. No.	Description
00-8909-007-00	Complete Cast Removal System

INCLUDES		
Product No.	Description	Packaging
00-8909-001-00	Zimmer Cast Cutter II	1
00-8909-066-00	Cast Cutter II Vacuum Connector Kit — Straight	1
00-8899-002-00	Cast Vacuum	1
00-8899-026-00	Zimmer Cast Cutter support frame	1

CAST VACUUM

- For removing cast dust and particles during cast cutting
- Helps keep cast room clean
- Wheel-mounted for easy transportation

Prod. No.	Description
00-8899-002-00	Cast Vacuum

CAST VACUUM ACCESSORIES		
Prod. No.	Description	Packaging
00-8909-015-00	Dust Collector/Blade Shield Assembly, 90-Degree	1
00-8909-016-00	Dust Collector/Blade Shield Assembly, Straight	1
00-8909-065-00	Vacuum Connector System Includes: 8909-015 90-Degree Shroud, 8899-518 Hose, 5625-380 Plastic Ties	
00-8909-066-00	Vacuum Connector System Includes: 8909-016 Straight, 8899-518 Hose, 5625-380 Plastic Ties	
00-8899-019-00	Vacuum Replacement Filter Bags	5 each
00-8899-026-00	Cast Cutter Support Frame	1



CASTING PRODUCTS

CAST CART

- Designed to hold materials for cast application and removal
- Enamel-coated steel cabinet with plastic work surface
- Side mounted cast bucket keeps water away from storage areas
- 6 buckets included



Prod. No.	Description	Dimensions
00-3820-000-00	Cast Cart	33 1/2 in. wide 35 1/2 in. high 18 1/2 in. deep (85cm x 47cm x 90cm)

CAST BUCKETS

- Wax covered double wrapped paper
- Disposable
- 5-quart capacity
- Sold in boxes of 25



Prod. No.	Description
00-3820-054-00	Cast Buckets (25)

CAST INSTRUMENTS

- Assortment of instruments for trimming cast padding and/or opening casts
- Stainless steel



00-3835-000-00



00-3807-002-00



00-3812-000-00



00-1720-000-00



00-3809-000-00



00-3826-000-00

Prod. No.	Description	Material	Dimensions
00-3835-000-00	Boehler Cast Breaker	Stainless steel	10 1/8 in. long (26cm) 3 1/4 in. jaw length (8mm) 2 1/2 in. max. jaw opening (6cm)
00-3807-002-00	Universal Scissors	Stainless steel with plastic handle	7 in. long (18cm)
00-3812-000-00	Brun Shears	Stainless steel	9 1/2 in. long (24cm) 1 1/4 in. sharp blade (3cm) 2 1/8 in. serrated blade (5cm)
00-1720-000-00	Esmarch Shears	Stainless steel	8 in. long (20cm) 2 in. blade length (5cm) 1/8 in. protective tip (3mm)
00-3809-000-00	Utility Felt Shears	Stainless steel	8 1/4 in. long (21cm) 2 3/8 in. & 2 11/16 in. blade lengths (6cm & 7cm)
00-3826-000-00	Cast Spreader	Stainless steel	8 3/4 in. long (22cm) 1 in. max. spread (3cm)

PADS AND POSITIONERS

ABDUCTION PILLOW—VINYL

- Maintains hip in abduction following surgery
- Vinyl-coated surface for easy cleaning
- Contact closure straps for easy application and adjustment



Prod. No.	Size	Proximal End Width	Distal End Width
00-1086-002-00	Small	7 3/16 in. (18cm)	14 in. (36cm)
00-1086-003-00	Medium	9 1/4 in. (23cm)	17 1/4 in. (44cm)
00-1086-004-00	Large	11 3/8 in. (29cm)	29 in. (74cm)

VINYL ABDUCTION PILLOW BANDS

- Replacement bands for 1086 abduction pillows
- Polyester elastic and cotton flannel



Prod. No.	Pillow	Size
00-1086-012-00	00-1086-002-00	Small
00-1086-013-00	00-1086-003-00	Medium
00-1086-014-00	00-1086-004-00	Large

VINYL ABDUCTION PILLOW STRAPS

- Replacement straps for 1086 abduction pillows
- Polyester pile webbing with contact closures



Prod. No.	Pillow	Size
00-1086-032-00	00-1086-002-00	Small
00-1086-033-00	00-1086-003-00	Medium
00-1086-034-00	00-1086-004-00	Large

ABDUCTION PILLOW TUBULAR PADS

- Replacement tubular pads for 1086 abduction pillows
- Packaged 4 per bag



Prod. No.	Pillow	Size
00-1086-022-00	00-1086-002-00	Small
00-1086-023-00	00-1086-003-00	Medium
00-1086-024-00	00-1086-004-00	Large



PADS AND POSITIONERS

ZIMFOAM® ABDUCTION PILLOW

- Supports and aligns patient's legs
- Contoured sides help support leg when patient is turned
- Long foam straps accommodate most calf or thigh sizes
- Disposable
- Sold in cases of 6
- Vacuum packed



Prod. No.	Depth	Length	Width
00-2671-001-00	6 in. (15cm)	17 1/4 in. (44cm)	12 in. (30cm)
00-2671-002-00	6 in. (15cm)	23 in. (58cm)	14 in. (36cm)
00-2671-003-00	6 in. (15cm)	26 in. (66cm)	16 in. (41cm)

BODY ALIGNER

- Grooved surface helps disburse pressure and enhances air circulation
- Supports side-lying patients at a 45 or 90 degree angle
- Sold in cases of 6
- Vacuum packed



Prod. No.	Dimensions
00-2665-000-00	10 in. x 10 in. x 22 in. (25cm x 25cm x 56cm)

CAST ELEVATOR

- For patients with leg casts, edema or those requiring leg elevation and rotation control
- Tunneled area can accommodate ice bags
- Supports leg in flexed or extended position
- Sold in cases of 6
- Vacuum packed



Prod. No.	Dimensions
00-2666-000-00	10 in. x 16 in. x 11 in. (25cm x 41cm x 28cm)

CRADLE BOOT

- Protects heels against pressure sores
- Supports foot in upright, neutral position
- Sold in cases of 6
- Vacuum packed



Prod. No.	Dimensions
00-2670-000-00	8 in. side (20cm) 14 in. long (36cm) 9 in. high (23cm)

PADS AND POSITIONERS

DEROTATION BOOT

- Hold foot firmly upright to guard against foot drop
- Straps hold boot securely to leg



Prod. No.	Description	Height	Length
00-2680-000-00	Complete Boot	11 1/4 in. (28cm)	8 1/2 in. (22cm)
00-2680-011-00	Liner		

FOOT DROP STOP

- Protects against pressure sores of the heel and ankle
- Helps hold foot in upright position
- Fully adjustable and disposable
- Packaged 6 per box



Prod. No.	Height	Length	Width
00-2669-000-00	6 in. (15cm)	12 in. (30cm)	8 in. (20cm)

LAMINECTOMY ARM CRADLE

- Supports patient's arm during surgery to help protect ulnar nerve
- Disposal
- Packaged 12 per box



Prod. No.	Length	Width
00-2663-000-00	24 in. (61cm)	15 in. (38cm)

MULTI-CRADLE POSITIONER

- Intended for prone head support during surgery
- Can also provide support for arm, leg, ankle or foot
- Packaged 12 per box



Prod. No.	Diameter	Height
00-2658-000-00	10 in. (25cm)	7 in. (18cm)

OR DONUT

- All purpose pressure dispersion device
- Packaged 12 per box



Prod. No.	Hole Diameter	Overall Diameter	Thickness
00-2656-000-00	4 in. (10cm)	8 in. (20cm)	3 in. (8cm)



PADS AND POSITIONERS

OR HEAD POSITIONER

- Scalloped foam construction helps disperse pressure on head during surgery
- Especially helpful during long surgical procedures
- Packaged 12 per box



Prod. No.	Height	Length	Width
00-2679-000-00	4 in. (10cm)	8 1/2 in. (22 cm)	8 in. (20cm)

ULNAR NERVE PROTECTOR

- Helps protect ulnar nerve during surgery
- Easy to apply with contact closure
- Packaged 12 per box



Prod. No.	Length
00-2664-000-00	8 in. (20cm)

L-P ULNAR NERVE PROTECTOR

- Holds arm off bed to help protect ulnar nerve
- Quickly applied with contact closure
- Economically priced
- Packaged 20 per box



Prod. No.	Length	Width
00-2678-000-00	6 in. (15cm)	8 in. (20cm)

ULTRA-CARE® HEEL/ELBOW PROTECTOR

- Helps reduce pressure to protect heels and elbows against pressure sores
- Tubular stockinette keeps sleeve in place with minimal constriction
- Foam pad molded to stockinette. No stitch seams to cause irritation
- Packaged 2 per box
- Washable



Prod. No.	Length
00-1782-001-00	9 in. (23cm)

PADS AND POSITIONERS

SLEEVE KNIT HEEL/ELBOW PROTECTOR

- Economical protection for heels and ankles
- Washable
- Packaged 10 per box



Prod. No.	Length
00-2790-010-00	7 1/4 in. (18cm)

ECONOMY HEEL/ELBOW PROTECTORS

- An economical means of decubitus protection
- Easily applied with contact closure
- Packaged 2 per box



Prod. No.	Length
00-1761-001-00	4 in. x 4 in. x 4 in. (10cm x 10cm x 10cm)

ZIMMER® HEEL/ELBOW PROTECTORS

- Polyester pile lining provides soft cushion
- Easily applied with contact closure
- Washable
- Packaged 2 per box



Prod. No.	Length	Width
00-0658-001-00	4 in. (10cm)	4 in. (10cm)



TRACTION BOOTS

UNIVERSAL TRACTION BOOT

- Universal size for inventory convenience
- Medial and lateral stays may be positioned anywhere on boot
- Excess material may be trimmed for proper fit



Prod. No.	Length	Max Calf Circumference
00-1774-001-00	20 in. (51cm)	20 in. (51cm)

DELUXE ZIM-TRAC® CONVOLUTED TRACTION BOOT

- Convolutions reduce ankle pressure and increase air circulation.
- Supports foot to help prevent foot drop
- Spreader block included



Prod. No.	Size	Length	Max Calf Circumference
00-2753-010-00	XS	13 in. (33cm)	6 5/8 in. - 11 3/4 in. (17cm - 30cm)
00-2753-011-00	S	20 in. (51cm)	8 3/8 in. - 13 3/4 in. (22cm - 35cm)
00-2753-012-00	M	20 in. (51cm)	9 5/8 in. - 15 3/4 in. (24cm - 40cm)
00-2753-013-00	L	20 in. (51cm)	10 7/8 in. - 17 3/4 in. (28cm - 45cm)
00-2753-014-00	XL	20 in. (51cm)	12 1/8 in. - 19 3/4 in. (31cm - 50cm)

ECONOMY LEG TRACTION BOOT

- Polyurethane convoluted foam with medial/lateral stays
- Does not include footplate or spreader block




Prod. No.	Size	Length	Max Calf Circumference
00-8712-012-00	S	20 in. (51cm)	8 3/8 in. - 13 3/4 in. (22cm - 35cm)
00-8712-013-00	M	20 in. (51cm)	9 5/8 in. - 15 3/4 in. (24cm - 40cm)
00-8712-014-00	L	20 in. (51cm)	10 7/8 in. - 17 3/4 in. (28cm - 45cm)
00-8712-015-00	XL	20 in. (51cm)	12 1/8 in. - 19 3/4 in. (31cm - 50cm)

PELVIC BELTS

ZIMCODE® PELVIC TRACTION BELTS

- Contact closure strap
- Wide range of sizes
- Belt width 8 1/2 in. (22cm)
- Double strap



Prod. No.	Size
00-0743-001-00	26 in. - 30 in. (66cm - 76cm)
00-0743-002-00	30 in. - 34 in. (76cm - 86cm)
00-0743-003-00	34 in. - 38 in. (86cm - 97cm)
00-0743-004-00	38 in. - 42 in. (97cm - 107cm)
00-0743-005-00	42 in. - 46 in. (107cm - 117cm)
00-0743-006-00	46 in. - 50 in. (117cm - 127cm)

ZIMMER® COMFORT PANEL PELVIC BELT

- Special elastic panels distribute support evenly for patient comfort
- Soft inner lining reduces skin irritation
- Belt width 9 in. (23cm)



Prod. No.	Size
00-0584-021-00	26 in. - 45 in. (66cm - 114cm)
00-0584-022-00	30 in. - 52 in. (76cm - 132cm)
00-0584-023-00	36 in. - 62 in. (91cm - 157cm)

UNIVERSAL SINGLE STRAP PELVIC TRACTION BELT


- Universal size for inventory convenience
- Accommodates wide range of hip sizes
- Contact closure
- Belt width 8 in. (20cm)



Prod. No.	Size
00-1678-001-00	26 in. - 48 in. (66cm - 122cm)

UNIVERSAL DOUBLE STRAP PELVIC BELT

- Straps with contact closure attachments can be positioned anywhere on belt to achieve most effective line of pull
- Universal size for inventory convenience
- Belt width 8 in. (20cm)



Prod. No.	Size
00-1570-001-00	26 in. - 48 in. (66cm - 122cm)
00-1570-002-00	48 in. - 60 in. (122cm - 152cm)

DELUXE FOAM PADDED TRACTION BELT

- Posterior strap directs maximum traction to pelvic region
- Contact closures for easy application and adjustment
- Belt width 8 1/2 in. (22cm)



Prod. No.	Size
00-0738-001-00	26 in. - 34 in. (66cm - 86cm)
00-0738-002-00	34 in. - 42 in. (86cm - 107cm)
00-0738-003-00	42 in. - 50 in. (107cm - 127cm)
00-0738-006-00	50 in. - 58 in. (127cm - 147cm)



HEAD HALTERS

NO. 7 HEAD HALTER

- Designed for intermittent overhead traction
- No spreader bar needed
- Ideal for clinic or home use
- Universal size
- Fully adjustable



Prod. No.

00-2700-007-00

ZIMMER® PRO HEAD HALTER

- Specifically designed for bed traction
- Directs pull to occipital region
- No spreader bar needed
- Universal size



Prod. No.

00-2700-003-00

DELUXE DISKARD HEAD HALTER

- Universal size
- Spreader bar required
- Foam padded for comfort
- Sold in boxes of 12



Prod. No.

00-0912-000-00

Packaging

12 per box

DELUXE PADDED HEAD HALTER

- Foam padded for patient comfort
- Contact closure for easy application
- Spreader bar included



Prod. No.

00-1778-001-00

Size

Small

00-1778-002-00

Medium

00-1778-003-00

Large



Please refer to package insert for complete product information, including contraindications, warnings, precautions, and adverse effects.

Contact your Zimmer Representative or visit us at www.zimmer.com.