



Zimmer®  
Pulsavac®  
Plus Wound  
Debridement  
Family



Superior options for maximum control





## Pulsavac® Plus Family

---

The *Pulsavac*® Plus, *Pulsavac* Plus LP and *Pulsavac* Plus AC wound debridement systems offer versatile solutions to wound irrigation, wound cleaning and therapy. Three models, featuring lightweight design and trigger-controlled pressure settings, offer a powerful cleaning action for bone or a gentle lavage for soft tissue. Concurrent suction and irrigation help to thoroughly remove debris and contaminants without flooding the field. The *Pulsavac* Plus and *Pulsavac* Plus AC are systems for use in orthopaedic and trauma cases. The *Pulsavac* Plus LP provides gentle lavage for soft tissue debridement.



**Tip Retaining Mechanism**

*provides a positive lock for secure attachment of tips.*

**Textured Surface**

*near tip to help provide a sure grip and firm control, even in awkward positions.*

**Rocker Switch**

*provides convenient pressure control at your fingertips.*

**Textured Pistol-Grip Handle**

*provides a secure, comfortable grip to facilitate control.*

**Ergonomic, Lightweight Design**

*fits comfortably into the palm of your hand facilitating use in tight areas.*

**Independent Fluid and Suction Tubes**

*provide concurrent fluid delivery and suction.*

**Power In Your Hands**

The Pulsavac Plus and Pulsavac Plus AC systems are ideal for use in the operating room and orthopaedic surgery applications. High power helps to effectively clean bone and cement debris during total joint procedures. Special tips were specifically designed for knee and hip procedures.



## THREE PULSAVAC PLUS MODELS AVAILABLE

### Pulsavac Plus •

Combines convenience and power for smooth, efficient pulsative irrigation and debridement. The external battery pack hangs outside the sterile field, and the variety of tips and trigger-controlled pressure offers either powerful cleaning or gentle lavage. The primary application site is the operating room for orthopaedic surgery. *Pulsavac Plus* operates on eight AA batteries.



**Reusable AC Adapter** eliminates the need for alkaline batteries, reducing the amount of waste.



### Pulsavac Plus AC •

An environmentally-friendly option. The reusable AC Adapter Power Pack must be purchased separately. This system completely eliminates the need for alkaline batteries.



### Pulsavac Plus LP •

This low pressure system delivers an adjustable, gentle lavage for wound cleaning and therapies and is an excellent alternative to traditional whirlpool therapy treatments. Primary application may occur in the OR, ER, wound care clinics, nursing homes or long-term care facilities. Therapy can be done in the patient's room or at the patient's home.

#### Low Pressure Debridement Tip Only



**00-5150-186-00** Flexible Tip, 10.5 inch length



## DEBRIDEMENT TIPS AVAILABLE

Debridement tips include the dual-cone splash shield, designed to minimize splashback and help control the dangers associated with translucent spray.

The color coded dots to the left of the tips indicate system compatibility.

● ● ● High Capacity Fan Spray w/Splash Shield, 12.7cm length  
 ● ● ● Without Splash Shield

● ● High Capacity Fan Spray w/Filter

● ● ● High Capacity Shower Spray w/Splash Shield, 12.7cm length  
 ● ● ● Without Splash Shield

● ● High Capacity Wide Single Stream, 7.2cm length

● ● High Capacity Four Orifices, 20.32cm x 1.27cm

● ● High Capacity Radial Spray, 22.86cm x .89cm

● ● High Capacity Brush Tip Radial Spray, 22.86cm x .89cm x 1.52cm

● ● High Capacity Hip Revision Tip, 23cm length

● ● High Capacity Four Orifices, 30.48cm x 1.27cm

● ● High Capacity Tibial Plateau Brush Tip, 7cm length

● ● High Capacity Adjustable Length w/Splash Shield

● ● High Capacity Single Stream, w/Splash Shield, 12.7cm length

# Ordering Information

## Pulsavac Plus Kits (Box 10)

Prod. No.	Description
00-5150-420-00	Pulsavac Plus Component Kit
00-5150-475-00	Pulsavac Plus Fan Kit
00-5150-482-00	Pulsavac Plus Hip Kit
00-5150-495-00	Pulsavac Plus Shower Kit

## Pulsavac Plus Kits (Box 1)

Prod. No.	Description
00-5150-420-01	Pulsavac Plus Component Kit
00-5150-475-01	Pulsavac Plus Fan Kit
00-5150-482-01	Pulsavac Plus Hip Kit
00-5150-495-01	Pulsavac Plus Shower Kit

## Pulsavac Plus AC Kits (Box 1)

Prod. No.	Description
00-5150-350-00	Pulsavac Plus AC Power Supply (US)
00-5150-426-01	Pulsavac Plus AC Component Kit
00-5150-476-01	Pulsavac Plus AC Fan Spray Kit
00-5150-486-01	Pulsavac Plus AC Hip Kit
00-5150-496-01	Pulsavac Plus AC Shower Spray Kit

## Pulsavac Plus LP Kits (Box 10)

Prod. No.	Description
00-5150-920-00	Pulsavac Plus LP Component Kit
00-5150-975-00	Pulsavac Plus LP Fan Spray Kit
00-5150-995-00	Pulsavac Plus LP Shower Spray Kit

## Multi-Orifice Tips (Box 10)

Prod. No.	Description
00-5150-175-00	High Capacity Fan Spray w/Splash Shield, 12.7cm length
00-5150-176-00	High Capacity Fan Spray w/Filter
00-5150-195-00	High Capacity Shower Spray w/Splash Shield, 12.7cm length
00-5150-375-00	High Capacity Fan Spray w/out Splash Shield, 10cm length
00-5150-395-00	High Capacity Shower Spray w/out Splash Shield, 10cm length

## Debridement/Suction Tips (Box 10)

Prod. No.	Description
00-5150-155-00	High Capacity Single Stream Tip w/Splash Shield, 12.7cm length
00-5150-065-00	Large Fluid Orifice Tip w/Wide Single Stream, 7.2cm length

## Femoral Tips (Box 10)

Prod. No.	Description
00-5150-181-00	High Capacity IM Tip Four Orifices, 20.32cm x 1.27cm
00-5150-182-00	High Capacity IM Tip w/Radial Spray, 22.86cm x .89cm
00-5150-183-00	High Capacity IM Brush Tip w/Radial Spray, 22.86cm x .89cm x 1.52cm
00-5150-184-00	High Capacity IM Hip Revision Tip, 23cm length
00-5150-185-00	High Capacity IM Tip Four Orifices, 30.48cm x 1.27cm
00-5150-187-00	High Capacity Tibial Plateau Brush Tip, 7cm length

## Specialty Tip (Box 10)

Prod. No.	Description
00-5150-145-00	High Capacity Adjustable Length Tip w/Splash Shield
00-5150-186-00	Flexible Tip, 10.5 inch length

## Miscellaneous (Box 10)

Prod. No.	Description
00-5150-188-00	18cm Splash Shield
00-5150-189-00	23cm Splash Shield

## AC Power Pack (Reusable) (Box 1)



Contact your Zimmer representative or visit us at [www.zimmer.com](http://www.zimmer.com)



+H124975150303001/\$081111L084